

## PRX II antibody, C-term

Cat. No. GTX89331

|              |            |
|--------------|------------|
| Host         | Goat       |
| Clonality    | Polyclonal |
| Isotype      | IgG        |
| Applications | WB         |
| Reactivity   | Human      |

Package  
100 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution                | Recommended dilution |
|-----------------------------------|----------------------|
| WB                                | 0.3-1µg/ml           |
| Not tested in other applications. |                      |

Calculated MW 22 kDa. ( [Note](#) )

Product Note This antibody is expected to recognise isoforms a and b only (NP\_005800.3 and NP\_859427.1 respectively).

## Properties


|               |                                                                                                                                                                                                                                                                                                                                      |
|---------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                                                                                                                               |
| Buffer        | TBS, 0.5% BSA                                                                                                                                                                                                                                                                                                                        |
| Preservative  | 0.02% Sodium azide                                                                                                                                                                                                                                                                                                                   |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.                                                                                                           |
| Concentration | 0.50 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                                                                                                                          |
| Immunogen     | Peptide with sequence C-NVDDSKKEYFSKHN, from the C Terminus of the protein sequence according to NP_005800.3; NP_859427.1.                                                                                                                                                                                                           |
| Purification  | Purified by ammonium sulphate precipitation followed by antigen affinity chromatography                                                                                                                                                                                                                                              |
| Conjugation   | Unconjugated                                                                                                                                                                                                                                                                                                                         |
| Note          | For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.<br><br>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product. |



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 06 Page 1 of 2

## DATA IMAGES



### GTX89331 WB Image

WB analysis of human brain lysate using GTX89331 PRX II antibody, C-term.

Dilution : 1µg/ml

Loading : 35µg protein in RIPA buffer



For full product information, images and publications, please visit our [website](#).