

SERPINI2 antibody, Internal

Cat. No. GTX89529

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.3-1µg/ml
IHC-P	3µg/ml

Note : Human Pancreas shows cytoplasmic staining of acinar cells and of rare cells in the islets of Langerhans.

Not tested in other applications.

Calculated MW	46 kDa. (Note)
Product Note	This antibody is expected to recognise both reported isoforms (NP_006208.1; NP_001012303.1).

Properties

Form	Liquid
Buffer	TBS, 0.5% BSA
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.50 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Peptide with sequence C-KGDWKQKFRKEDTQ, from the internal region of the protein sequence according to NP_006208.1; NP_001012303.1.
Purification	Purified by ammonium sulphate precipitation followed by antigen affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTx89529 WB Image

WB analysis of human pancreas lysate using GTx89529 SERPINI2 antibody, Internal.

Dilution : 0.3µg/ml

Loading : 35µg protein in RIPA buffer


GTx89529 IHC-P Image

IHC-P analysis of human pancreas using GTx89529 SERPINI2 antibody, Internal.

Antigen retrieval : citrate buffer pH 6

Dilution : 3µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).