

FGFR1 antibody, C-term

Cat. No. GTX89588

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.3-1µg/ml
IHC-P	4µg/ml

Note : Human Small Intestine shows staining of the enterocytic cytoplasm.

Not tested in other applications.

Calculated MW	92 kDa. (Note)
Product Note	This antibody is expected to recognise five isoforms (as represented by NP_075598.2; NP_056934.2; NP_075593.1; NP_075594.1; NP_075599.1).

Properties

Form	Liquid
Buffer	TBS, 0.5% BSA
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.50 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Peptide with sequence CLPRHPAQLANGGLKR, from the C Terminus of the protein sequence according to NP_075598.2; NP_056934.2; NP_075593.1; NP_075594.1; NP_075599.1.
Purification	Purified by ammonium sulphate precipitation followed by antigen affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTX89588 IHC-P Image

IHC-P analysis of human small intestine using GTX89588 FGFR1 antibody, C-term.

Antigen retrieval : citrate buffer pH 6

Dilution : 4µg/ml




GTX89588 WB Image

WB analysis of human breast lysate using GTX89588 FGFR1 antibody, C-term.

Dilution : 0.3µg/ml

Loading : 35µg protein in RIPA buffer



GTX89588 IHC-P Image

IHC-P analysis of human kidney using GTX89588 FGFR1 antibody, C-term.

Antigen retrieval : citrate buffer pH 6

Dilution : 4µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).