

Triosephosphate isomerase antibody, C-term

Cat. No. GTX89594

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.01-0.03µg/ml
IHC-P	5µg/ml

Note : Human Liver shows staining of the sinusoids.

Not tested in other applications.

Calculated MW	31 kDa. (Note)
Product Note	This antibody is expected to recognise reported isoforms 1 and 2 (NP_000356.1; NP_001152759.1).

Properties

Form	Liquid
Buffer	TBS, 0.5% BSA
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.50 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Peptide with sequence C-LKPEFVDIINAKQ, from the C Terminus of the protein sequence according to NP_000356.1; NP_001152759.1.
Purification	Purified by ammonium sulphate precipitation followed by antigen affinity chromatography
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 08 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTx89594 IHC-P Image

IHC-P analysis of human pancreas using GTX89594 Triosephosphate isomerase antibody, C-term.
Antigen retrieval : Tris/EDTA buffer pH 9
Dilution : 5µg/ml



GTx89594 WB Image

WB analysis of human liver lysate using GTX89594 Triosephosphate isomerase antibody, C-term.
Dilution : 0.003µg/ml
Loading : 30µg protein in RIPA buffer



GTx89594 IHC-P Image

IHC-P analysis of human liver using GTX89594 Triosephosphate isomerase antibody, C-term.
Antigen retrieval : Tris/EDTA buffer pH 9
Dilution : 2µg/ml



For full product information, images and publications, please visit our [website](#).