

LXR alpha / LXR beta antibody, Internal

Cat. No. GTX89656

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Mouse

References (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1-3µg/ml

Not tested in other applications.

Product Note

This antibody is expected to recognise an epitope corresponding to aa 429-442 of human LXR alpha protein (NP_005684.2) and to aa 443-456 of human LXR beta protein (NP_009052.3).

Properties

Form	Liquid
Buffer	TBS, 0.5% BSA
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.50 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Peptide with sequence CRLQDKKLPLLSEI, from the internal region of the protein sequence according to NP_005684.2; NP_009052.3.
Purification	Purified by ammonium sulphate precipitation followed by antigen affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTx89656 WB Image

WB analysis of human spleen lysate using GTx89656 LXR alpha / LXR beta antibody, Internal.


Dilution : 1µg/ml

Loading : 35µg protein in RIPA buffer



GTx89656 WB Image

The data was published in the Front Vet Sci in 2020. [PMID: 33195535](https://pubmed.ncbi.nlm.nih.gov/33195535/)



GTx89656 WB Image

The data was published in the Front Vet Sci in 2020. [PMID: 33195535](https://pubmed.ncbi.nlm.nih.gov/33195535/)



For full product information, images and publications, please visit our [website](https://www.genetex.com).