

Retinoid X Receptor beta antibody, Internal

Cat. No. GTX89670

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human

References (1)

Package

100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1-3µg/ml

Not tested in other applications.

Calculated MW 57 kDa. ([Note](#))

Product Note

This antibody is expected to recognise an epitope corresponding to aa 70-83 of human RXR beta protein and does not cross-react with either RXR alpha or gamma.

Properties

Form	Liquid
Buffer	TBS, 0.5% BSA
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.50 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Peptide with sequence CRDGMGDSGRDSRSP, from the internal region of the protein sequence according to NP_068811.1.
Purification	Purified by ammonium sulphate precipitation followed by antigen affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 08 Page 1 of 2

DATA IMAGES

GTX89670 WB Image

WB analysis of A431 lysate using GTX89670 Retinoid X Receptor beta antibody, Internal.

Dilution : 1µg/ml

Loading : 35µg protein in RIPA buffer


GTX89670 WB Image

WB analysis of HEK293 overexpressing Retinoid X Receptor beta (mock transfection in first lane) using GTX89670 Retinoid X Receptor beta antibody, Internal.



For full product information, images and publications, please visit our [website](https://www.genetex.com).