

DLC1 (Isoforms 1 and 3) antibody, N-term

Cat. No. GTX89689

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	5µg/ml

Note : Human Adrenal Gland (Medulla) shows strong reticulate staining in the cytoplasm of medular cells.

Not tested in other applications.

Calculated MW	171 kDa. (Note)
Product Note	This antibody is expected to recognise isoforms 1 and 3 (NP_872584.2 and NP_079043.3 resp) but NOT isoform 2 (NP_006085.2).

Properties

Form	Liquid
Buffer	TBS, 0.5% BSA
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.50 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Peptide with sequence SVAIRKRSWEEHC, from the N Terminus of the protein sequence according to NP_872584.2; NP_079043.3.
Purification	Purified by ammonium sulphate precipitation followed by antigen affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX89689 IHC-P Image

IHC-P analysis of human adrenal gland using GTX89689 DLC1 (Isoforms 1 and 3) antibody, N-term.

Antigen retrieval : citrate buffer pH 6

Dilution : 5µg/ml



For full product information, images and publications, please visit our [website](#).