

RNF36 antibody, C-term

Cat. No. GTX89832

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.5-2µg/ml
IHC-P	2.5-3.8µg/ml

Note : Human Testis shows speckled staining in nuclei of spermatocytes and spermatids.

Not tested in other applications.

Calculated MW	57 kDa. (Note)
Product Note	This antibody is expected to recognize both reported isoforms (NP_892030.3; NP_542783.2).

Properties

Form	Liquid
Buffer	TBS, 0.5% BSA
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.50 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Peptide with sequence C-GENKEPLHILHPQ, from the C Terminus of the protein sequence according to NP_892030.3; NP_542783.2.
Purification	Purified by ammonium sulphate precipitation followed by antigen affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 06 Page 1 of 2

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTx89832 WB Image

WB analysis of mouse testis lysate using GTx89832 RNF36 antibody, C-term.

Dilution : 0.5µg/ml

Loading : 35µg protein in RIPA buffer


GTx89832 IHC-P Image

IHC-P analysis of human testis using GTx89832 RNF36 antibody, C-term.

Antigen retrieval : citrate buffer pH 6

Dilution : 3.8µg/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).