

SOCS4 blocking peptide

Cat. No. GTX89892-PEP

Species	Human
----------------	-------

Package	100 µg
----------------	--------

Applications

Product Note	The peptide was used in the production of GTX89892.
---------------------	-----------------------------------------------------

Properties

Form	Lyophilized powder
-------------	--------------------

Buffer	Reconstitute with 200µl distilled water to obtain a 0.5mg/ml peptide solution. Lyophilized from 5% (v/v) acetonitrile/H ₂ O.
---------------	-----------------------------------------------------------------------------------------------------------------------------------------

Preservative	No preservatives
---------------------	------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
----------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Region/Sequence	C-SKVRVLRIDAPEQQ
------------------------	------------------

Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
-------------	---------------------------------------------------------------------------------------------------------------------------------------------------

Note	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.
-------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------



For full product information, images and publications, please visit our [website](#).