

DAP Kinase 2 antibody, C-term

Cat. No. GTX89970

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1-3µg/ml
IHC-P	2-4µg/ml

Note : Human Breast shows staining of the luminal side of the cell membranes of epithelial cells in the ducts.

Not tested in other applications.

Calculated MW 43 kDa. ([Note](#))

Properties


Form	Liquid
Buffer	TBS, 0.5% BSA
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.50 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Peptide with sequence C-KALHPRRRSSTS, from the C Terminus of the protein sequence according to NP_055141.2.
Purification	Purified by ammonium sulphate precipitation followed by antigen affinity chromatography
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 07 Page 1 of 2

DATA IMAGES




GTX89970 IHC-P Image

IHC-P analysis of human breast using GTX89970 DAP Kinase 2 antibody, C-term.

Antigen retrieval : citrate buffer pH 6

Dilution : 2.5µg/ml




GTX89970 WB Image

WB analysis of mouse brain extracts using GTX89970 DAP Kinase 2 antibody, C-term.

Dilution : 2µg/ml

Loading : 35µg protein in RIPA buffer



GTX89970 IHC-P Image

IHC-P analysis of paraffin embedded human thyroid using GTX89970 DAP Kinase 2 antibody, C-term.

Antigen retrieval : Steamed antigen retrieval with citrate buffer Ph 6

Dilution : 2.5µg/ml



For full product information, images and publications, please visit our [website](#).