

## K-Ras antibody [4F3]

Cat. No. GTX91809

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, ELISA
Reactivity	Human, Cat

## Package

100 µg

## Applications

## Application Note

For ELISA (tested against immunogen): Use at an assay dependent dilution. For WB (tested against immunogen): Use at an assay dependent dilution. Optimal dilutions/concentrations should be determined by the researcher.

**Calculated MW** 22 kDa. ( [Note](#) )**Product Note** This antibody may cross react with N-Ras and H-Ras proteins.

## Properties

**Form** Liquid**Buffer** Phosphate-buffered saline, pH 7.2, containing no preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.5 mg/ml (Please refer to the vial label for the specific concentration.)**Immunogen** KRAS (NP\_004976, 16 a.a. ~ 126 a.a) partial recombinant protein**Purification** Protein A affinity purified**Conjugation** Unconjugated

**Note**

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).