



Accelerate Your Oncology Research Breakthroughs

Your Expertise

Our Antibodies

Accelerated Discovery



GeneTex's High-Quality
Recombinant Antibodies
for Oncology Research

Why Choose GeneTex Recombinant Antibodies for Your Research?



Reproducible

Eliminating batch-to-batch variation ensures reliable data across studies.



High Sensitivity and Specificity

Engineered for precise target recognition, delivering exceptionally accurate data.




Rigorous Validation

Follows stringent "5 + 1 Pillar" validation protocols, with an emphasis on KO/KD characterization.




Carrier-Free Buffer Formulation


Ready for conjugation and other advanced applications.




Claudin 18.2 antibody [HL5230]
(GTX645486)



B-Raf (V600E Mutant) antibody
[HL5673] (GTX646482)




CEACAM1 antibody [HL5159]
(GTX645355)




L1CAM antibody [HL5786]
(GTX646857)

Cat. No.	Product Name	Host & Clonality	Reactivity	Applications
GTX637428	14-3-3 antibody [HL1776]	Rec Rb mAb	Hu	WB
GTX646219	14-3-3 zeta antibody [HL5545]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX640930	4E-BP1 (phospho Thr37/Thr46) antibody [HL3276]	Rec Rb mAb	Hu, Ms	WB
GTX635479	53BP1 antibody [HL275]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX638580	Actin antibody [HL2372]	Rec Rb mAb	AGMK, Cat, Dm, Dog, Hu, Ms, Rat, Rb	WB, ICC/IF, IHC-P
GTX636668	Activin Receptor Type 1B antibody [HL1266]	Rec Rb mAb	Hu, Ms	WB
GTX648176	ACVRL1 antibody [HL6091]	Rec Rb mAb	Ms, Rat	WB
GTX643635	ADAM10 antibody [HL4378]	Rec Rb mAb	Hu	IHC-P
GTX641694	ADAR1 antibody [HL3666]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P
GTX638758	Adenosine A1 Receptor antibody [HL2442]	Rec Rb mAb	Dog, Hu, Ms, Rat	WB, IHC-P
GTX642209	Adenosine A2b Receptor antibody [HL3959]	Rec Rb mAb	Ms	IHC-Fr
GTX638125	ADFP antibody [HL2147]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF
GTX648333	ADNP antibody [HL6241]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX637225	AF9 antibody [HL1644]	Rec Rb mAb	Hu	WB, ICC/IF
GTX638885	AGTR1 antibody [HL2524]	Rec Rb mAb	Hu, Ms	WB, ICC/IF, IHC-P
GTX637885	AHR antibody [HL1984]	Rec Rb mAb	Hu, Ms	WB, ICC/IF
GTX640148	AKT (phospho Ser473) antibody [HL2871]	Rec Rb mAb	Hu, Ms	WB, IHC-P
GTX641363	AKT (phospho Thr450) antibody [HL3478]	Rec Rb mAb	Hu, Ms	WB
GTX640258	AKT antibody [HL2915]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX640359	AKT1 (phospho Ser473) antibody [HL2972]	Rec Rb mAb	Hu, Ms	WB
GTX636413	AKT1 antibody [HL1142]	Rec Rb mAb	Cat, Dog, Hu, Rat	WB, ICC/IF, IHC-P
GTX648159	Aldolase A antibody [HL6080]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX648328	Aldolase B antibody [HL6236]	Rec Rb mAb	Ms, Rat	WB, IHC-P
GTX644904	ALIX antibody [HL4841]	Rec Rb mAb	Hu, Ms	WB, ICC/IF
GTX636670	Alkaline phosphatase (intestinal) antibody [HL1268]	Rec Rb mAb	Dm, Hu, Ms	WB, IHC-P
GTX637904	Alkaline Phosphatase (Tissue Non-Specific) antibody [HL2003]	Rec Rb mAb	Dog, Hu, Ms	WB, IHC-P
GTX637965	ALKBH5 antibody [HL2061]	Rec Rb mAb	Hu, Ms, Rat	WB
GTX637983	ALP / PDLIM3 antibody [HL2075]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX646675	alpha 1 Antitrypsin antibody [HL5724]	Rec Rb mAb	Hu	WB, IHC-P
GTX642060	alpha Actinin 4 antibody [GT95]	Rec Ms mAb	Hu	WB
GTX637291	alpha Actinin antibody [HL1692]	Rec Rb mAb	Hu, Ms	WB, ICC/IF, IHC-P
GTX648204	Alpha fetoprotein / AFP antibody [HL6119]	Rec Rb mAb	Hu	WB, ICC/IF
GTX646520	alpha Internexin antibody [HL5689]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX638140	alpha Tubulin antibody [HL2162]	Rec Rb mAb	AGMK, Cat, Dog, E. coli, ... ¹	WB, ICC/IF, IHC-P
GTX642037	ALY antibody [HL3813]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX640142	AMFR antibody [HL2865]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX642027	AMIGO1 antibody [GT94]	Rec Ms mAb	Ms, Rat	WB
GTX636639	Androgen Receptor (ARv7 Splice Variant) antibody [HL1240]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P, IHC-P (cell pellet)
GTX636021	Androgen Receptor antibody [HL1049]	Rec Rb mAb	Hu, Rat	WB, ICC/IF, IHC-P
GTX638218	Angiogenin antibody [HL2214]	Rec Rb mAb	Hu, Ms	WB, ICC/IF, IHC-P
GTX642156	Angiopoietin 2 antibody [HL3911]	Rec Rb mAb	Hu	IHC-P
GTX647569	Anillin antibody [HL5957]	Rec Rb mAb	Hu	WB, ICC/IF
GTX638119	Annexin A1 antibody [HL2141]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX644038	APC antibody [HL4473]	Rec Rb mAb	Hu	WB, ICC/IF
GTX637312	Apolipoprotein A1 antibody [HL1713]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, ELISA
GTX648181	Aquaporin 4 antibody [HL6096]	Rec Rb mAb	Ms	WB, IHC-P
GTX640697	ARL13B antibody [HL3179]	Rec Rb mAb	Hu, Ms	WB, ICC/IF
GTX644466	ARL5A antibody [HL4602]	Rec Rb mAb	Hu	ICC/IF, IHC-P
GTX638542	Aromatase antibody [HL2339]	Rec Rb mAb	Hu	WB, IHC-P
GTX636575	Artemis antibody [HL1227]	Rec Rb mAb	Hu	WB, ICC/IF


¹ AGMK, Cat, Dog, E. coli, Hm, Hu, Ms, Rat, Rb, Zfsh




RAS (G12D Mutant) antibody [HL10] (GTX635362)



14-3-3 zeta antibody [HL5545] (GTX646219)






Her2 / ErbB2 (phospho Tyr1139) antibody [HL3272] (GTX640926)




EGFRvIII antibody [HL2541] (GTX638910)


Cat. No.	Product Name	Host & Clonality		Reactivity	Applications
GTX641873	ASK1 (phospho Ser966) antibody [HL3701]	Rec Rb	mAb	Hu	WB, ICC/IF
GTX641552	ATF1 (phospho Ser63) antibody [HL3581]	Rec Rb	mAb	Hu	WB, ICC/IF
GTX641410	ATF4 antibody [HL3513]	Rec Rb	mAb	Ms	WB
GTX638551	ATG4B antibody [HL2348]	Rec Rb	mAb	Hu, Ms, Rat	WB, IHC-P
GTX648822	ATIC antibody [HL6383]	Rec Rb	mAb	Hu, Ms	WB
GTX636086	ATM (phospho Ser1981) antibody [HL1062]	Rec Rb	mAb	Hu	WB, ICC/IF, IHC-P
GTX635395	ATM antibody [2C1-RB]	Rec Rb	mAb	Hu, Ms, Rat	WB, IHC-P
GTX642199	ATP synthase gamma antibody [HL3949]	Rec Rb	mAb	Hu, Ms, Rat	WB, ICC/IF
GTX643629	ATP5D antibody [HL4372]	Rec Rb	mAb	Hu, Ms, Rat, Zfsh	WB, IHC-P
GTX640318	ATP5J antibody [HL2934]	Rec Rb	mAb	Hu, Ms, Rat	WB, IHC-P
GTX636734	ATP6V1H antibody [HL1316]	Rec Rb	mAb	Dm, Hu, Ms	WB, IHC-P
GTX642327	ATP7A antibody [GT128]	Rec Ms	mAb	Hu	WB
GTX640889	ATR (phospho Ser428) antibody [HL3246]	Rec Rb	mAb	Hu, Ms	WB, ICC/IF
GTX635571	ATR (phospho Thr1989) antibody [HL132]	Rec Rb	mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX635478	ATRX antibody [HL1210]	Rec Rb	mAb	Hu	WB, ICC/IF
GTX640260	Aurora A antibody [HL2917]	Rec Rb	mAb	Hu	WB, ICC/IF
GTX644810	Aurora B antibody [HL4798]	Rec Rb	mAb	Hu	WB, ICC/IF
GTX647608	Axin 2 antibody [HL5994]	Rec Rb	mAb	Hu	WB
GTX648397	AXIN1 antibody [HL6295]	Rec Rb	mAb	Hu	WB
GTX642023	AXL antibody [HL3804]	Rec Rb	mAb	Hu, Ms, Rat	WB
GTX643004	BAI1 antibody [HL4190]	Rec Rb	mAb	Hu, Ms, Rat	WB, IHC-P
GTX642532	BAI2 antibody [HL4093]	Rec Rb	mAb	Hu, Ms, Rat	WB
GTX642534	BAI3 antibody [HL4095]	Rec Rb	mAb	Hu, Ms, Rat	WB, ICC/IF
GTX637619	BATF antibody [HL1870]	Rec Rb	mAb	Hu, Ms	WB, ICC/IF
GTX635715	Bax antibody [HL236]	Rec Rb	mAb	Chinese Hamster, Dog, Hu, Ms, Rat	WB, IHC-P
GTX638176	BCA1 antibody [HL2172]	Rec Rb	mAb	Hu	WB, IHC-P
GTX639471	BCKDK antibody [HL2710]	Rec Rb	mAb	Hu, Ms, Rat	WB, IHC-P
GTX647192	BCL10 antibody [HL5868]	Rec Rb	mAb	Hu, Ms, Rat	WB
GTX638196	BCL6 antibody [HL2192]	Rec Rb	mAb	Dog, Hu, Ms	WB
GTX637940	Bcl-XL antibody [HL2039]	Rec Rb	mAb	Cat, Dog, Hu	WB, ICC/IF, IHC-P
GTX645327	BDKRB2 antibody [HL5131]	Rec Rb	mAb	Hu	WB
GTX637675	beta Actin antibody [HL1926]	Rec Rb	mAb	AGMK, Cat, Dm, Dog, Hu, Ms, Rat, Rb, Zfsh	WB, ICC/IF, IHC-P
GTX645299	beta Catenin (phospho Ser33) antibody [HL5105]	Rec Rb	mAb	Hu, Ms, Rat	WB
GTX641356	beta Catenin (phospho Ser675) antibody [HL3471]	Rec Rb	mAb	Hu, Ms	WB
GTX633010	beta Catenin antibody [GT2169]	Rec Ms	mAb	Dog, Hu, Ms	WB, ICC/IF, IHC-P
GTX645319	beta Tubulin antibody [HL5123]	Rec Rb	mAb	Hu, Ms	WB, ICC/IF
GTX647571	BIN3 antibody [HL5959]	Rec Rb	mAb	Hu, Ms, Rat	WB
GTX639636	BLM antibody [HL2771]	Rec Rb	mAb	Hu	WB
GTX648170	BMP4 antibody [HL6085]	Rec Rb	mAb	Hu	WB, ICC/IF
GTX636515	BNIP3L antibody [HL1203]	Rec Rb	mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX646156	BOB1 antibody [HL5492]	Rec Rb	mAb	Hu, Ms	WB
GTX642339	B-Raf (G469A Mutant) antibody [HL4008]	Rec Rb	mAb	Hu	WB, ICC/IF, IHC-P (cell pellet)
GTX642341	B-Raf (L597R Mutant) antibody [HL4010]	Rec Rb	mAb	Hu	WB, IHC-P (cell pellet)
GTX646482	B-Raf (V600E Mutant) antibody [HL5673]	Rec Rb	mAb	Hu	IHC-P
GTX637085	B-Raf antibody [HL1606]	Rec Rb	mAb	Cat, Dog, Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX636696	BRD9 antibody [HL1282]	Rec Rb	mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX642128	BTK (phospho Tyr223) antibody [HL3887]	Rec Rb	mAb	Hu, Ms	WB
GTX646464	BTK (phospho Tyr551) antibody [HL5655]	Rec Rb	mAb	Hu	WB
GTX637455	BTK antibody [HL1803]	Rec Rb	mAb	Hu, Ms	WB, ICC/IF
GTX641407	C Reactive Protein antibody [GT55]	Rec Ms	mAb	Hu, Ms, Rat	WB, ELISA, Sandwich ELISA




Cat. No.	Product Name	Host & Clonality	Reactivity	Applications
GTX636508	  C1 inhibitor antibody [HL1196]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX640102	   C3a Receptor antibody [HL2831]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P, FCM
GTX644789	  CADM4 antibody [GT272]	Rec Ms mAb	Rat	WB
GTX638641	  Calceineurin B antibody [HL2419]	Rec Rb mAb	Hu, Ms, Rat	WB
GTX646257	  Caldesmon antibody [HL5575]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF
GTX637077	   Calnexin antibody [HL1598]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	WB, IHC-P
GTX638748	  Calpain 5 antibody [HL2432]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX637961	   Calreticulin antibody [HL2057]	Rec Rb mAb	Hu	WB
GTX641653	  CaMKII (phospho Thr286) antibody [HL3662]	Rec Rb mAb	Hu	WB
GTX636876	  Carbonic Anhydrase IX antibody [HL1410]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX638312	   CARD10 antibody [HL2260]	Rec Rb mAb	Hu, Ms	WB, IHC-P
GTX638136	  CASC5 antibody [HL2158]	Rec Rb mAb	Hu	WB, ICC/IF
GTX637906	    Caspase 1 antibody [HL2005]	Rec Rb mAb	Hu	WB
GTX636810	    Caspase 3 antibody [HL1370]	Rec Rb mAb	Hu, Rat	WB, IHC-P
GTX640257	    Caspase 8 antibody [HL2914]	Rec Rb mAb	Hu	WB
GTX642096	  Caspase 9 (cleaved Asp330) antibody [HL3858]	Rec Rb mAb	Hu	WB
GTX642095	  Caspase 9 antibody [HL3857]	Rec Rb mAb	Hu	WB
GTX648319	  Catalase antibody [HL6227]	Rec Rb mAb	Hu, Rat	WB
GTX638369	    Cathepsin S antibody [HL2302]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX648113	  Caveolin 1 antibody [HL6039]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX639554	    CCN3 antibody [HL2726]	Rec Rb mAb	Hu, Rat	WB, ICC/IF, IHC-P
GTX640472	    CCR1 antibody [HL3034]	Rec Rb mAb	Cat, Hu	WB, ICC/IF, IHC-P
GTX640109	  CCR2 antibody [HL2838]	Rec Rb mAb	Hu	ICC/IF
GTX638830	    CCR4 antibody [HL2481]	Rec Rb mAb	Hu, Ms, Rat	ICC/IF, IHC-P, FCM
GTX640520	  CCR7 antibody [HL3074]	Rec Rb mAb	Hu	ICC/IF, IHC-P
GTX640497	  CCR9 antibody [HL3054]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P, FCM
GTX640123	   CCR10 antibody [HL2852]	Rec Rb mAb	Hu	ICC/IF, Functional Assay
GTX640848	  CCR11 antibody [HL3208]	Rec Rb mAb	Ms	IHC-P
GTX643044	   CCRL2 antibody [HL4222]	Rec Rb mAb	Hu	WB
GTX647183	  CD2 antibody [HL5859]	Rec Rb mAb	Hu	WB, IHC-P
GTX645318	  CD2AP antibody [HL5122]	Rec Rb mAb	Hu, Ms, Rat	WB
GTX647567	  CD3 delta antibody [HL5955]	Rec Rb mAb	Hu, Ms	WB, ICC/IF, IHC-P
GTX647166	  CD3 epsilon antibody [HL5848]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX648193	  CD5 antibody [HL6108]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX645504	  CD7 antibody [HL5248]	Rec Rb mAb	Hu	WB, IHC-P
GTX644095	   CD9 antibody [HL4525]	Rec Rb mAb	Hu, Rat	WB, IHC-P, ELISA, Sandwich ELISA
GTX640543	  CD11b antibody [HL3097]	Rec Rb mAb	Ms, Rat	WB, IHC-P
GTX648344	  CD13 antibody [HL6252]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX645784	  CD19 antibody [HL5319]	Rec Rb mAb	Hu	WB, ICC/IF
GTX648311	  CD20 antibody [HL6219]	Rec Rb mAb	Hu	WB, IHC-P
GTX645787	  CD22 antibody [HL5322]	Rec Rb mAb	Hu	WB
GTX643014	   CD23 antibody [GT176]	Rec Ms mAb	Hu	WB, ICC/IF
GTX637068	  CD28 antibody [HL1589]	Rec Rb mAb	Cat, Dog, Hu	WB, ICC/IF, IHC-P
GTX645788	  CD30 antibody [HL5323]	Rec Rb mAb	Hu	WB, ICC/IF
GTX642317	    CD31 antibody [HL3991]	Rec Rb mAb	Hu, Ms	WB, ICC/IF, IHC-P, IHC-Fr
GTX637785	    CD32B antibody [HL1948]	Rec Rb mAb	Hu, Rat	WB, ICC/IF, IHC-P
GTX645077	   CD33 antibody [HL4991]	Rec Rb mAb	Hu	WB, ICC/IF
GTX646688	  CD38 antibody [HL5737]	Rec Rb mAb	Hu	WB, ICC/IF
GTX648933	  CD40 antibody [HL6485]	Rec Rb mAb	Hu	WB
GTX637956	  CD41 antibody [HL2052]	Rec Rb mAb	Hu	WB




CD155 antibody [HL5788]
(GTX646859)



FGL1 antibody [HL2317]
(GTX638480)




CD28 antibody [HL1589]
(GTX637068)




Galectin 3 antibody [HL2396]
(GTX638613)


Cat. No.	Product Name	Host & Clonality	Reactivity	Applications
GTX645276	CD43 antibody [HL5086]	Rec Rb mAb	Hu, Ms	WB, ICC/IF, IHC-P
GTX637231	CD44 antibody [HL1650]	Rec Rb mAb	Cat, Dog, Hu	WB, ICC/IF, IHC-P, FCM, IP
GTX648359	CD47 antibody [HL6267]	Rec Rb mAb	Hu, Ms	WB, ICC/IF, IHC-P
GTX643611	CD63 antibody [HL4361]	Rec Rb mAb	Hu, Ms, Cat	WB, IHC-P, ELISA
GTX638816	CD66e antibody [HL2467]	Rec Rb mAb	Hu	WB, IHC-P
GTX648089	CD70 antibody [HL6015]	Rec Rb mAb	Hu	WB
GTX646846	CD71 antibody [HL5775]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX638851	CD73 antibody [HL2493]	Rec Rb mAb	Hu, Ms, Rat, Zfsh	WB, IHC-P
GTX641627	CD74 antibody [HL3636]	Rec Rb mAb	Hu	WB, ICC/IF
GTX646282	CD79a antibody [HL5600]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P
GTX637265	CD81 antibody [HL1667]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P, FCM, ELISA
GTX637416	CD90 antibody [HL1766]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX646689	CD99 antibody [HL5738]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P
GTX648454	CD105 antibody [HL6344]	Rec Rb mAb	Hu	WB, IHC-P
GTX637124	CD133 antibody [HL1631]	Rec Rb mAb	Hu	FCM
GTX646859	CD155 antibody [HL5788]	Rec Rb mAb	Hu	WB, IHC-P
GTX642518	CD169 antibody [HL4079]	Rec Rb mAb	Hu, Ms	ICC/IF, IHC-P
GTX637310	CD252 antibody [HL1711]	Rec Rb mAb	Hu	WB
GTX638307	CDC25C antibody [HL2255]	Rec Rb mAb	Hu	WB
GTX635978	CDH11 antibody [HL1044]	Rec Rb mAb	Cat, Hu, Ms, Rat	WB
GTX641639	CDK1 / CDC2 (phospho Thr161) + CDK2 (phospho Thr160) + CDK3 (phospho Thr160) antibody [HL3648]	Rec Rb mAb	Hu	WB
GTX636948	CDK2 (phospho Tyr15) antibody [HL1474]	Rec Rb mAb	Dm, Hu, Ms, Rat	WB, IHC-P
GTX636873	CDK4 antibody [HL1407]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX637986	CDK6 antibody [HL2078]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX645355	CEACAM1 antibody [HL5159]	Rec Rb mAb	Hu	WB, IHC-P
GTX644085	CELSR1 antibody [HL4515]	Rec Rb mAb	Hu	WB
GTX639628	CEP152 antibody [HL2763]	Rec Rb mAb	Hu	WB, ICC/IF
GTX637970	CHD2 antibody [HL2066]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX637450	CHD4 antibody [HL1798]	Rec Rb mAb	Hu, Ms	WB, IHC-P
GTX639462	CHI3L1 antibody [HL2708]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P
GTX643560	Chk1 (phospho Ser317) antibody [HL4313]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF
GTX635572	Chk1 (phospho Ser345) antibody [HL122]	Rec Rb mAb	Hu	WB, ICC/IF
GTX635979	cIAP1 antibody [HL1045]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	WB, IHC-P
GTX639474	c-Jun (phospho Ser73) antibody [HL2713]	Rec Rb mAb	Hu	WB, ICC/IF
GTX639112	c-Kit (D816V Mutant) antibody [HL2653]	Rec Rb mAb	Hu	WB, ICC/IF
GTX645342	c-Kit antibody [HL5146]	Rec Rb mAb	Hu, Rat	WB, ICC/IF
GTX646212	Claudin 1 antibody [HL5538]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX637000	Claudin 3 antibody [HL1527]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P
GTX640666	Claudin 5 antibody [HL3156]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX645116	Claudin 18 antibody [HL5023]	Rec Rb mAb	Hu, Ms	WB
GTX645486	Claudin 18.2 antibody [HL5230]	Rec Rb mAb	Hu, Ms	WB, IHC-P, IHC-P (cell pellet)
GTX640244	c-Met (phospho Tyr1234) antibody [HL2907]	Rec Rb mAb	Hu	WB
GTX647222	c-Met antibody [HL5898]	Rec Rb mAb	Hu	WB, ICC/IF
GTX647191	C-Myb antibody [HL5867]	Rec Rb mAb	Hu, Ms	WB
GTX648180	c-Myc antibody [HL6095]	Rec Rb mAb	Hu	WB, ICC/IF
GTX639037	CNOT8 antibody [HL2586]	Rec Rb mAb	Hu, Japanese Medaka, Ms, Rat, Zfsh	WB, ICC/IF
GTX641584	Cofilin 1 antibody [HL3609]	Rec Rb mAb	Hu, Ms	WB, ICC/IF
GTX640483	COL1A1 antibody [HL3040]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF
GTX637953	COL1A2 antibody [HL2049]	Rec Rb mAb	Hu, Ms, Rat	WB




HIF1 alpha antibody [HL3154]
(GTX640664)



LDHA antibody [HL2615]
(GTX639069)




VEGFA antibody [HL1755]
(GTX637405)




eNOS (phospho Ser1177)
antibody [HL3498] (GTX641386)

Cat. No.	Product Name	Host & Clonality	Reactivity	Applications
GTX636771	COL4A1 antibody [HL1351]	Rec Rb mAb	Hu, Ms	WB, ICC/IF, IHC-P
GTX638633	Collagen I + Collagen II + Collagen III antibody [HL2048 + HL1907]	Rec Rb MultiAb	Hu, Ms, Rat	WB, IHC-P
GTX637424	Collagen II antibody [HL1772]	Rec Rb mAb	Hu	WB
GTX637655	Collagen III antibody [HL1906]	Rec Rb mAb	Hu, Ms	WB, ICC/IF
GTX640855	COMP antibody [HL3215]	Rec Rb mAb	Hu	ELISA
GTX638316	COX4 antibody [HL2264]	Rec Rb mAb	Hu, Ms, Rat, Zfsh	WB, ICC/IF, IHC-P
GTX638475	CP110 antibody [HL2312]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX642305	CPN1 antibody [HL3979]	Rec Rb mAb	Hu, Zfsh	WB
GTX636468	CPT1A antibody [HL1164]	Rec Rb mAb	Cat, Dog, Hu, Ms	WB, IHC-P
GTX639641	CREB (phospho Ser133) antibody [HL2776]	Rec Rb mAb	Hu, Ms	WB, IHC-P
GTX640158	CREB antibody [HL2879]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX642114	CRLR antibody [HL3876]	Rec Rb mAb	Hu	WB
GTX639560	CTLA4 antibody [HL2732]	Rec Rb mAb	Hu	WB, IHC-P
GTX644702	CX3CR1 antibody [HL4701]	Rec Rb mAb	Hu	ICC/IF, FCM
GTX639338	CXCR1 antibody [HL2674]	Rec Rb mAb	Hu, Ms, Rat	IHC-P
GTX639056	CXCR2 antibody [HL2604]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P, FCM, Functional Assay
GTX639064	CXCR4 antibody [HL2612]	Rec Rb mAb	Dog, Hu, Ms, Rat	WB, ICC/IF, IHC-P, FCM
GTX641880	CXCR5 antibody [HL3708]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P, IHC-Fr
GTX640120	CXCR6 antibody [HL2849]	Rec Rb mAb	Hu	WB, ICC/IF
GTX638193	CXCR7 antibody [HL2189]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P, IHC-Fr
GTX642109	CXCL10 / IP10 antibody [HL3871]	Rec Rb mAb	Hu	IHC-P
GTX632502	CXCL16 antibody [GT516]	Rec Ms mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX638223	Cyclin D1 antibody [HL2219]	Rec Rb mAb	Dog, Hu, Rat	WB
GTX645275	CYFIP2 antibody [HL5085]	Rec Rb mAb	Hu, Ms	WB, ICC/IF
GTX636662	CYLD antibody [HL1260]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX638122	CYR61 antibody [HL2144]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX638209	Cytochrome C antibody [HL2205]	Rec Rb mAb	Cat, Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX648153	Cytokeratin 7 antibody [HL6074]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P
GTX636469	Cytokeratin 19 antibody [HL1165]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX641071	Cytokeratin 20 antibody [HL3309]	Rec Rb mAb	Hu, Ms	WB, IHC-P
GTX644955	DARC antibody [HL4890]	Rec Rb mAb	Hu	ICC/IF
GTX641696	Desmin antibody [HL3668]	Rec Rb mAb	Hu, Ms	WB, IHC-P
GTX644508	Desmocollin 2 antibody [HL4655]	Rec Rb mAb	Hu	WB, ICC/IF
GTX637937	Desmoglein 2 antibody [HL2036]	Rec Rb mAb	Hu, Ms	WB, IHC-P
GTX636580	Dicer antibody [HL1232]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX70141	DNA ligase I antibody [10H5]	Rec Ms mAb	Hu, Ms	WB, ICC/IF, IHC-P, IP
GTX638333	DNA ligase III antibody [HL2280]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX636986	DNA polymerase beta antibody [HL1512]	Rec Rb mAb	Hu, Japanese Medaka, Ms, Rat, Zfsh	WB, ICC/IF, IHC-P
GTX643622	DNA polymerase gamma antibody [HL4365]	Rec Rb mAb	Hu, Ms	WB, ICC/IF
GTX641072	DNA polymerase iota antibody [HL3310]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX643030	DNA-PKcs (phospho Ser2056) antibody [HL4207]	Rec Rb mAb	Hu	WB, ICC/IF
GTX638541	DUSP8 antibody [HL2338]	Rec Rb mAb	Hu, Rat	WB, ICC/IF
GTX636527	DUT antibody [HL1210]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P
GTX648101	E2F3 antibody [HL6027]	Rec Rb mAb	Hu	WB, ICC/IF
GTX640841	E-Cadherin (phospho Ser838/Ser840) antibody [HL3201]	Rec Rb mAb	Hu, Ms, Rat	WB
GTX636577	E-Cadherin antibody [HL1229]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX641652	EDG4 / LPAR2 antibody [HL3661]	Rec Rb mAb	Hu	ICC/IF
GTX635716	EGFR (E746-A750del) antibody [HL363]	Rec Rb mAb	Hu	WB
GTX635834	EGFR (L858R Mutant) antibody [HL1022]	Rec Rb mAb	Hu	WB
GTX645769	EGFR (phospho Ser1047) antibody [HL5304]	Rec Rb mAb	Hu	WB




CD31 antibody [HL3991]
(GTX642317)



CD105 antibody [HL6344]
(GTX648454)




PDGF Receptor alpha antibody
[HL5001] (GTX645094)



Thrombospondin 2 antibody
[HL1140] (GTX636411)

Cat. No.	Product Name	Host & Clonality	Reactivity	Applications
GTX640912	EGFR (phospho Thr669) antibody [HL3265]	Rec Rb mAb	Hu	WB, ICC/IF
GTX641131	EGFR (phospho Tyr1045) antibody [HL3354]	Rec Rb mAb	Hu	WB, ICC/IF
GTX640938	EGFR (phospho Tyr1086) antibody [HL3284]	Rec Rb mAb	Hu	WB, ICC/IF
GTX640932	EGFR (phospho Tyr1148) antibody [HL3278]	Rec Rb mAb	Hu	WB, ICC/IF
GTX640934	EGFR (phospho Tyr1173) antibody [HL3280]	Rec Rb mAb	Hu	WB, ICC/IF
GTX641096	EGFR (phospho Tyr845) antibody [HL3331]	Rec Rb mAb	Hu	WB, ICC/IF
GTX641365	EGFR (phospho Tyr974) antibody [HL3480]	Rec Rb mAb	Hu	WB
GTX642543	EGFR antibody [HL1097]	Rec Rb mAb	Hu, Ms	WB, ICC/IF, IHC-P
GTX638910	EGFRvIII antibody [HL2541]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P (cell pellet)
GTX645438	EGR1 antibody [HL5186]	Rec Rb mAb	Hu	WB, ICC/IF
GTX637028	eIF4E antibody [HL1553]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX637408	EIF6 antibody [HL1758]	Rec Rb mAb	Dm, Hu, Japanese Medaka, Ms, Rat, Zfsh	WB, IHC-P
GTX643709	EMR2 antibody [HL4420]	Rec Rb mAb	Hu	WB, ICC/IF
GTX648195	EMR3 antibody [HL6110]	Rec Rb mAb	Hu	IHC-P
GTX641386	eNOS (phospho Ser1177) antibody [HL3498]	Rec Rb mAb	Hu	WB
GTX636759	EpCAM antibody [HL1339]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P, IHC-Fr, FCM
GTX639110	EphA3 antibody [HL2651]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX637446	EPO antibody [HL1794]	Rec Rb mAb	Cat, Dog, Hu	WB, IHC-P, ELISA, ... ¹
GTX635617	EERK1 (phospho Thr202/Tyr204) + ERK2 (phospho Thr185/Tyr187) antibody [HL173]	Rec Rb mAb	Dm, Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX640331	ERK1/2 antibody [HL2946]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX646181	ESAM antibody [HL5518]	Rec Rb mAb	Hu	WB, ICC/IF
GTX636202	Estrogen Receptor alpha antibody [HL1081]	Rec Rb mAb	Dog, Hu	WB, ICC/IF, IHC-P
GTX638311	ETS2 antibody [HL2259]	Rec Rb mAb	Hu, Ms	WB, ICC/IF
GTX641453	Exonuclease 1 antibody [HL3547]	Rec Rb mAb	Hu	WB, ICC/IF
GTX640252	Ezrin (phospho Thr567) + Radixin (Thr564) + Moesin (Thr558) Antibody [HL2909]	Rec Rb mAb	Hu	WB
GTX638490	Ezrin antibody [HL2327]	Rec Rb mAb	Dog, Hu, Ms, Rat	WB
GTX646110	F4/80 antibody [HL5468]	Rec Rb mAb	Hu, Ms	WB, ICC/IF, IHC-P, FCM
GTX637863	FAH antibody [HL1970]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX643064	FAK (phospho Tyr407) antibody [HL4231]	Rec Rb mAb	Hu, Ms	WB
GTX635714	FAK (phospho Tyr576) antibody [HL127]	Rec Rb mAb	Hu, Ms, Rat	WB
GTX648928	FAK antibody [HL6480]	Rec Rb mAb	Hu, Ms	WB
GTX646839	Fascin 1 antibody [HL5768]	Rec Rb mAb	Hu, Ms	WB, ICC/IF
GTX636884	FCER1G antibody [HL1418]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	WB, ICC/IF, IHC-P, IP
GTX648318	FEN1 antibody [HL6226]	Rec Rb mAb	Ms	WB, IHC-P
GTX641372	FGFR1 (phospho Tyr766) + FGFR2 (phospho Tyr769) + FGFR3 (phospho Tyr760) antibody [HL3487]	Rec Rb mAb	Hu, Ms	WB
GTX641587	FGFR3 (phospho Tyr724) antibody [HL3612]	Rec Rb mAb	Hu	WB
GTX637074	FGFR3 antibody [HL1595]	Rec Rb mAb	Cat, Dog, Hu	WB
GTX638480	FGL1 antibody [HL2317]	Rec Rb mAb	Hu	WB, ICC/IF
GTX648392	Fibrillin 1 antibody [HL6290]	Rec Rb mAb	Hu	WB
GTX639626	Fibronectin antibody [HL2761]	Rec Rb mAb	Hu	WB, ICC/IF
GTX637457	FLI1 antibody [HL1805]	Rec Rb mAb	Hu	WB, IP
GTX640353	Flightless I antibody [HL2966]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX643486	Flotillin 1 antibody [HL4250]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX648115	FLRT3 antibody [GT396]	Rec Ms mAb	Hu, Ms	WB, ICC/IF
GTX640686	FOXA1 antibody [HL3168]	Rec Rb mAb	Hu, Ms, Zfsh	WB, ICC/IF, IHC-P, IHC-Wm
GTX640131	FOXO1 (phospho Thr24) antibody [HL2860]	Rec Rb mAb	Hu	WB
GTX640166	FOXO1 (phospho Ser256) antibody [HL2887]	Rec Rb mAb	Hu	IHC-P


¹ WB, IHC-P, ELISA, Neutralizing /Inhibition, Sandwich ELISA




E-Cadherin antibody [HL1229]
(GTX636577)



MMP9 antibody [HL3239]
(GTX640882)




MUC2 antibody [HL1724]
(GTX637372)



ZEB1 antibody [HL2245]
(GTX638294)


Cat. No.	Product Name	Host & Clonality	Reactivity	Applications
GTX638194	FOXO3A antibody [HL2190]	Rec Rb mAb	Hu, Ms	WB
GTX644721	FOXP3 antibody [HL4717]	Rec Rb mAb	Cat, Hu, Ms	WB, IHC-P
GTX643498	FPR1 antibody [HL4262]	Rec Rb mAb	Hu	ICC/IF, IHC-P, FCM
GTX644700	Frizzled 1 antibody [HL4699]	Rec Rb mAb	Hu	WB
GTX646694	Frizzled 2 antibody [HL5743]	Rec Rb mAb	Hu, Ms	WB
GTX640095	Frizzled 4 antibody [HL2824]	Rec Rb mAb	Hu, Zfsh	WB
GTX644920	Frizzled 6 antibody [HL4857]	Rec Rb mAb	Hu	WB
GTX648189	Frizzled 7 antibody [HL6104]	Rec Rb mAb	Hu	WB
GTX643662	Frizzled 8 antibody [HL4405]	Rec Rb mAb	Ms	WB
GTX637274	Frizzled 9 antibody [HL1675]	Rec Rb mAb	Hu	WB
GTX638772	FUS antibody [HL2454]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX648457	FUT4 antibody [HL6347]	Rec Rb mAb	Hu	WB
GTX645351	FXR1 antibody [HL5155]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF
GTX638303	G3BP1 antibody [HL2254]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX638314	GADD153 antibody [HL2262]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P
GTX638613	Galectin 3 antibody [HL2396]	Rec Rb mAb	Dog, Hu, Ms, Rat	WB, IHC-P, IHC-Fr, ELISA, Sandwich ELISA
GTX638601	GALNT6 antibody [HL2392]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX636480	gamma Tubulin antibody [HL1175]	Rec Rb mAb	AGMK, Cat, Dm, Dog, Hm, Hu, Ms, Rat, Rb, Zfsh	WB, ICC/IF, IHC-Fr
GTX636896	Gasdermin D antibody [HL1430]	Rec Rb mAb	Hu	WB, IHC-P
GTX648858	GATA6 antibody [HL6413]	Rec Rb mAb	Hu	WB
GTX641067	G-CSF antibody [GT21]	Rec Ms mAb	Hu	WB, ELISA, Sandwich ELISA
GTX643528	GDF15 antibody [HL4288]	Rec Rb mAb	Hu	IHC-P
GTX641355	GEF-H1 antibody [HL3470]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P, IP
GTX637767	Gelsolin antibody [HL1930]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF
GTX641902	Gephyrin antibody [HL3712]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX639685	GERP antibody [HL2797]	Rec Rb mAb	Cat, Hu, Rat	WB
GTX644234	GHS-R1 antibody [HL4566]	Rec Rb mAb	Hu, Ms	WB, IHC-P
GTX640943	GTR antibody [HL3289]	Rec Rb mAb	Hu	WB
GTX635619	Gli1 antibody [HL247]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P, IP
GTX643519	Glypican 1 antibody [HL4279]	Rec Rb mAb	Hu	WB, IHC-P
GTX645443	GM130 antibody [HL5191]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P
GTX642007	Golgin 97 antibody [HL3788]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF
GTX648257	GPNMB antibody [HL6170]	Rec Rb mAb	Hu	WB
GTX637589	GPRC5A antibody [HL1864]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P
GTX641566	GPRC5C antibody [HL3595]	Rec Rb mAb	Hu, Ms	IHC-P, IHC-Fr
GTX642471	GPRC5D antibody [HL4056]	Rec Rb mAb	Hu	WB, ICC/IF
GTX638329	GPS2 antibody [HL2276]	Rec Rb mAb	Hu, Ms, Rat, Zfsh	WB, ICC/IF, IHC-P
GTX647163	GPX1 antibody [HL5845]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF
GTX635477	GPX2 antibody [HL125]	Rec Rb mAb	Dog, Hu, Ms, Rat, Zfsh	WB, IHC-P
GTX642192	Granzyme B antibody [GT111]	Rec Ms mAb	Hu	WB
GTX637540	Grp78 antibody [HL1815]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX636766	GSPT1 antibody [HL1346]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX642155	HAPLN2 antibody [GT107]	Rec Ms mAb	Ms	WB, IHC-P
GTX647152	HAX1 antibody [HL5834]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P
GTX641720	hCG beta antibody [GT75]	Rec Ms mAb	Hu	ELISA, Lateral Flow, Sandwich ELISA
GTX637272	HCK antibody [HL1673]	Rec Rb mAb	Hu, Ms	WB, IHC-P
GTX647137	HE4 antibody [HL5819]	Rec Rb mAb	Hu	WB
GTX640829	Hec1 (phospho Ser 55) antibody [HL3190]	Rec Rb mAb	Hu	WB
GTX643045	Hec1 antibody [HL4223]	Rec Rb mAb	Hu	WB
GTX645508	HELIOS antibody [HL5252]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P




LI-Cadherin antibody [HL5636]
(GTX646325)



Claudin 1 antibody [HL5538]
(GTX646212)




Occludin antibody [HL6239]
(GTX648331)




Vimentin antibody [HL1506]
(GTX636980)


Cat. No.	Product Name	Host & Clonality	Reactivity	Applications
GTX637432	Heme Oxygenase 1 antibody [HL1780]	Rec Rb mAb	Cat, Hu, Ms, Rat	WB, IHC-P
GTX641464	Her2 / ErbB2 (phospho Tyr877) antibody [HL3558]	Rec Rb mAb	Hu	WB, ICC/IF
GTX640926	Her2 / ErbB2 (phospho Tyr1139) antibody [HL3272]	Rec Rb mAb	Hu	WB, ICC/IF
GTX641093	Her2 / ErbB2 (phospho Tyr1221/Tyr1222) antibody [HL3328]	Rec Rb mAb	Hu	WB
GTX636940	Her2 / ErbB2 (phospho Tyr1248) antibody [HL1466]	Rec Rb mAb	Hu, Rat	WB, ICC/IF, IHC-P
GTX641472	Her2 / ErbB2 antibody [HL3566]	Rec Rb mAb	Hu	WB, FCM
GTX638573	HEY1 antibody [HL2367]	Rec Rb mAb	Dog, Hu, Ms, Rat, Zfsh	WB, IHC-P
GTX637378	HGS antibody [HL1730]	Rec Rb mAb	Cat, Dog, Hu	WB, ICC/IF, IHC-P
GTX640664	HIF1 alpha antibody [HL3154]	Rec Rb mAb	Dog, Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX641412	Histone H2A.X antibody [HL3515]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX636713	Histone H2A.XS139ph (phospho Ser139) antibody [HL1299]	Rec Rb mAb	Hu, Ms	WB, ICC/IF, IHC-P
GTX637631	HLA class I antibody [HL1882]	Rec Rb mAb	Hu	WB, IHC-P
GTX639675	HLA-A2 antibody [BB7.2]	Rec Ms mAb	Hu	FCM
GTX639674	HLA-B7 antibody [BB7.1]	Rec Ms mAb	Hu	FCM
GTX637773	HMGA2 antibody [HL1936]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX645278	hnRNP A1 antibody [HL5088]	Rec Rb mAb	Hu	WB
GTX648882	HPD antibody [HL6437]	Rec Rb mAb	Hu, Ms, Rat	WB
GTX637440	HSC70 antibody [HL1788]	Rec Rb mAb	Hu, Japanese Medaka, Ms, Rat, Zfsh	WB, IHC-P
GTX641960	HSP27 (phospho Ser82) antibody [HL3761]	Rec Rb mAb	Hu	WB, ICC/IF
GTX637059	Hsp70 antibody [HL1580]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX644012	ICAM1 / CD54 antibody [HL4447]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P
GTX646193	ICAM3 antibody [HL5530]	Rec Rb mAb	Hu	WB, ICC/IF
GTX646294	ICAM4 antibody [HL5612]	Rec Rb mAb	Hu	WB, ICC/IF
GTX645732	ICAM5 antibody [HL5294]	Rec Rb mAb	Ms, Rat	WB, IHC-P
GTX644949	IDH1 (R132H Mutant) antibody [HL4884]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P, IHC-P (cell pellet)
GTX638817	IFIT3 antibody [HL2468]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P
GTX636979	IGFBP3 antibody [HL1505]	Rec Rb mAb	Hu	WB
GTX637052	IKK gamma antibody [HL1577]	Rec Rb mAb	Hu, Ms	WB, ICC/IF, IHC-P
GTX645733	IKZF3 antibody [HL5295]	Rec Rb mAb	Hu	WB
GTX641947	IL1 alpha antibody [HL3748]	Rec Rb mAb	Hu	WB
GTX646488	IL1 Receptor 2 antibody [HL5679]	Rec Rb mAb	Hu	WB
GTX640234	IL10 antibody [HL2897]	Rec Rb mAb	Hu	WB, ELISA
GTX647605	IL2 Receptor alpha antibody [HL5991]	Rec Rb mAb	Hu	WB, ICC/IF
GTX646192	IL13 Receptor alpha 1 antibody [HL5529]	Rec Rb mAb	Hu	WB
GTX645975	IL13 Receptor alpha 2 antibody [HL5407]	Rec Rb mAb	Hu, Ms	WB
GTX645246	IL2 antibody [GT322]	Rec Ms mAb	Hu	ELISA, Sandwich ELISA
GTX641122	IL5 antibody [GT30]	Rec Ms mAb	Hu	ELISA
GTX638598	IL6 antibody [HL2389]	Rec Rb mAb	Hu	WB, ICC/IF, ELISA, Sandwich ELISA
GTX637584	IL18 antibody [HL1859]	Rec Rb mAb	Hu, Ms	WB, ICC/IF, IHC-P, ELISA, Sandwich ELISA
GTX644912	IL21 antibody [HL4849]	Rec Rb mAb	Hu	ELISA
GTX636531	iNOS antibody [HL1213]	Rec Rb mAb	Ms, Rat	WB, ICC/IF, IHC-P
GTX645346	INPPL1 antibody [HL5150]	Rec Rb mAb	Hu	WB
GTX636833	Insulin antibody [HL1393]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P, IHC-Fr, ELISA
GTX648341	Integrin alpha 4 antibody [HL6249]	Rec Rb mAb	Hu	WB, ICC/IF
GTX641573	Integrin alpha V antibody [HL3598]	Rec Rb mAb	Hu, Rat	WB, IHC-P
GTX641905	Integrin beta 1 / CD29 antibody [HL3715]	Rec Rb mAb	Cat, Hu, Ms, Rat	WB, ICC/IF, IHC-P, IHC-Fr
GTX642035	Interferon gamma Receptor 1 antibody [HL3811]	Rec Rb mAb	Hu	WB
GTX637015	IP3 Receptor I antibody [HL1541]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX638293	IPMK antibody [HL2244]	Rec Rb mAb	Hu, Ms, Rat	WB
GTX644210	IRE1 alpha (phospho Ser724) antibody [HL4542]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF




ATM (phospho Ser1981) antibody [HL1062] (GTX636086)   



Cyclin D1 antibody [HL2219] (GTX638223)   




p53 antibody [HL2242] (GTX638291)   




mTOR antibody [HL2216] (GTX638220)   


Cat. No.	Product Name	Host & Clonality	Reactivity	Applications
GTX648292	  IRF3 antibody [HL6200]	Rec Rb mAb	Hu	WB
GTX647596	  IRF4 / MUM1 antibody [HL5982]	Rec Rb mAb	Hu	WB, ICC/IF
GTX637918	   ISG15 antibody [HL2017]	Rec Rb mAb	Cat, Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX641549	  ISG20 antibody [HL3578]	Rec Rb mAb	Hu, Ms	WB
GTX642194	  ITGBL1 antibody [HL3944]	Rec Rb mAb	Hu	WB
GTX636666	  ITK antibody [HL1264]	Rec Rb mAb	Hu, Ms	WB, ICC/IF, IHC-P
GTX642976	 JIP2 antibody [GT170]	Rec Ms mAb	Hu, Ms, Rat	WB
GTX635799	   JNK (phospho Thr183/Tyr185) antibody [HL1008]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX648273	  KDM4A antibody [HL6181]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX642563	   Ki67 antibody [HL4117]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P
GTX640486	  Kiss1 antibody [HL3043]	Rec Rb mAb	Hu	WB, IHC-P
GTX638321	   KLF5 antibody [HL2269]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX637048	 KPNA3 antibody [HL1573]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX644460	 Kv1.5 antibody [GT241]	Rec Ms mAb	Ms, Rat	WB
GTX646857	  L1CAM antibody [HL5786]	Rec Rb mAb	Hu, Ms, Rat	WB
GTX642177	  LAG3 antibody [HL3932]	Rec Rb mAb	Hu	WB
GTX640865	  Lamin B1 antibody [HL3224]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX648972	  Laminin beta 3 antibody [HL6516]	Rec Rb mAb	Hu	WB
GTX648973	  Laminin gamma 2 antibody [HL6517]	Rec Rb mAb	Hu	WB
GTX638556	   LAT1 / SLC7A5 antibody [HL2353]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX639686	  LC3B antibody [HL2798]	Rec Rb mAb	Hu, Ms, Rat, Zfsh	WB, ICC/IF, IHC-P
GTX639069	  LDHA antibody [HL2615]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P
GTX646679	  LEF1 antibody [HL5728]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX644815	  Leukotriene B4 Receptor 2 antibody [HL4803]	Rec Rb mAb	Hu	WB
GTX638146	   LGALS3BP antibody [HL2168]	Rec Rb mAb	Hu, Ms	WB, ICC/IF, IHC-P
GTX646325	  LI-Cadherin antibody [HL5636]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P
GTX648833	  LONP1 antibody [HL6394]	Rec Rb mAb	Hu, Ms	WB
GTX642555	  LOXL2 antibody [HL4109]	Rec Rb mAb	Hu	WB
GTX641354	  MAG antibody [HL3469]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P, IHC-Fr
GTX638090	   MAP3K8 antibody [HL2118]	Rec Rb mAb	Hu, Ms	WB, IHC-P
GTX642520	  MATH2 antibody [HL4081]	Rec Rb mAb	Ms, Rat	IHC-P
GTX641548	  MAX antibody [HL3577]	Rec Rb mAb	Cat, Hu, Ms, Rat	WB
GTX640674	  Maxi Potassium channel alpha antibody [HL3164]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF
GTX636945	   MC1 Receptor antibody [HL1471]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P, IHC-Fr
GTX640637	  MC4 Receptor antibody [HL3145]	Rec Rb mAb	Hu, Ms	IHC-P, FCM
GTX640696	  MC5 Receptor antibody [HL3178]	Rec Rb mAb	Hu, Ms	WB
GTX647174	  MCAK antibody [HL5856]	Rec Rb mAb	Hu	WB
GTX645272	  MCH-1R antibody [HL5082]	Rec Rb mAb	Ms	IHC-P
GTX637019	  MCL1 antibody [HL1544]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	WB, IHC-P
GTX642511	  MCM7 antibody [HL4072]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF
GTX645423	  MDM2 antibody [HL5171]	Rec Rb mAb	Hu	WB, ICC/IF
GTX647559	  MDMX / MDM4 antibody [HL5947]	Rec Rb mAb	Hu, Ms	WB
GTX642469	  MEK1 (phospho Ser217) + MEK2 (phospho Ser221) antibody [HL4054]	Rec Rb mAb	Hu, Ms	WB
GTX644736	   Melatonin Receptor 1A antibody [HL4732]	Rec Rb mAb	Hu, Ms, Rat	WB
GTX636729	  MFF antibody [HL1311]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX640350	  MIF antibody [HL2963]	Rec Rb mAb	Hu, Ms	WB, ELISA
GTX637437	  Mig-6 antibody [HL1785]	Rec Rb mAb	Hu	WB, IHC-P
GTX647156	  MITF antibody [HL5838]	Rec Rb mAb	Hu	WB, IHC-P, IHC-P (cell pellet)
GTX645198	  MLH1 antibody [HL5028]	Rec Rb mAb	Hu	WB, ICC/IF
GTX638279	  MLK4 antibody [HL2234]	Rec Rb mAb	Hu	WB, ICC/IF




CD44 antibody [HL1650]
(GTX637231)    



EpCAM antibody [HL1339]
(GTX636759)    




CD133 antibody [HL1631]
(GTX637124)    




SOX2 antibody [HL1193]
(GTX636505)    


Cat. No.	Product Name	Host & Clonality	Reactivity	Applications
GTX637574	   MLKL antibody [HL1849]	Rec Rb mAb	Hu, Rat	WB, ICC/IF
GTX636661	 MMP1 antibody [HL1259]	Rec Rb mAb	Hu	WB, ICC/IF
GTX636524	   MMP2 antibody [HL1207]	Rec Rb mAb	Hu, Ms	WB, IHC-P
GTX640882	  MMP9 antibody [HL3239]	Rec Rb mAb	Hu	WB, IHC-P
GTX640513	  MMP13 antibody [HL3067]	Rec Rb mAb	Hu, Ms	WB
GTX636826	   Mre11 antibody [HL1386]	Rec Rb mAb	Dog, Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX644746	  MRGPRG antibody [HL4742]	Rec Rb mAb	Hu	WB
GTX643710	  MRGPRX4 antibody [HL4421]	Rec Rb mAb	Hu	WB
GTX647585	  MSH2 antibody [HL5971]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF
GTX638918	   MSH3 antibody [HL2546]	Rec Rb mAb	Hu	WB
GTX645260	  MSH6 antibody [HL5070]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P
GTX641585	  MST1 (phospho Thr183) + MST2 (phospho Thr180) antibody [HL3610]	Rec Rb mAb	Hu	WB
GTX643023	   MTHFD2 antibody [HL4200]	Rec Rb mAb	Hu	WB
GTX638220	   mTOR antibody [HL2216]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P, IP
GTX640514	  mtRNA polymerase antibody [HL3068]	Rec Rb mAb	Hu	WB
GTX637913	   mtTFA antibody [HL2012]	Rec Rb mAb	Hu	WB, IHC-P
GTX648869	  MUC1 antibody [HL6424]	Rec Rb mAb	Hu	WB
GTX637372	   MUC2 antibody [HL1724]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX636960	 Mucin 4 antibody [HL1486]	Rec Rb mAb	Hu, Ms	WB, IHC-P
GTX636922	  MURF1 antibody [HL1450]	Rec Rb mAb	Hu, Ms	WB
GTX637955	   MX1 antibody [HL2051]	Rec Rb mAb	Cat, Dog, Hu, Ms	WB, ICC/IF, IHC-P
GTX641576	  MyD88 antibody [HL3601]	Rec Rb mAb	Hu	WB
GTX641208	  Myeloperoxidase antibody [HL3388]	Rec Rb mAb	Hu	WB
GTX643642	  MYL7 antibody [HL4385]	Rec Rb mAb	Ms, Rat	WB
GTX636812	   MyoD1 antibody [HL1372]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX637662	   Myogenin antibody [HL1913]	Rec Rb mAb	Hu, Ms, Zfsh	WB
GTX640170	  Myosin Light Chain 2 (MLC-2v) (phospho Ser15) antibody [HL2891]	Rec Rb mAb	Hu, Ms	WB, ICC/IF
GTX648199	  Napsin A antibody [HL6114]	Rec Rb mAb	Hu	WB, IHC-P
GTX640612	  NDP52 antibody [HL3126]	Rec Rb mAb	Hu	WB, IHC-P
GTX640923	  NDUFS2 antibody [HL3269]	Rec Rb mAb	Hu, Ms	WB
GTX648124	  Nectin 2 antibody [HL6045]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P
GTX642949	   NEIL2 antibody [HL4144]	Rec Rb mAb	Hu	WB, IHC-P
GTX638839	  Nek1 antibody [HL2489]	Rec Rb mAb	Hu, Ms, Rat	WB
GTX638872	  NEK2 antibody [HL2511]	Rec Rb mAb	Hu	WB, IHC-P
GTX636768	   NEK7 antibody [HL1348]	Rec Rb mAb	Hu, Ms	WB, ICC/IF, IHC-P
GTX648968	  Nestin antibody [HL6512]	Rec Rb mAb	Hu	WB
GTX644461	 Neurofascin antibody [GT242]	Rec Ms mAb	Rat	WB
GTX643062	   NFkB p65 antibody [HL4229]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF
GTX637436	   NFkB p105 antibody [HL1784]	Rec Rb mAb	Hu	WB
GTX648856	  NIPA antibody [HL6411]	Rec Rb mAb	Ms, Rat	WB
GTX639545	  NKG2A antibody [HL2718]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P
GTX639954	   NLRP3 antibody [HL2813]	Rec Rb mAb	Hu, Ms	WB, IP
GTX641070	  NRF1 antibody [HL3308]	Rec Rb mAb	Hu, Ms, Rat, Zfsh	WB, ICC/IF
GTX638901	   Nucleobindin 1 antibody [HL2532]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX642089	  Nucleophosmin (phospho Thr237) antibody [HL3855]	Rec Rb mAb	Hu, Ms, Rat	WB
GTX647229	  Nucleophosmin antibody [HL5905]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX636570	  Nucleoporin p62 antibody [HL1226]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX639324	  NUR77 antibody [HL2660]	Rec Rb mAb	Hu	WB
GTX638616	   OASL antibody [HL2399]	Rec Rb mAb	Hu, Rat	WB, IHC-P
GTX648331	  Occludin antibody [HL6239]	Rec Rb mAb	Hu, Ms	WB




CD81 antibody [HL1667]
(GTX637265)   



TSG101 antibody [4A10]
(GTX70255)    




CD63 antibody [HL4361]
(GTX643611)   




CD9 antibody [HL4525]
(GTX644095)   

Cat. No.	Product Name	Host & Clonality	Reactivity	Applications
GTX644243	   OLFML3 antibody [HL4567]	Rec Rb mAb	Hu	WB, IHC-P
GTX641197	   OPG antibody [HL3377]	Rec Rb mAb	Hu	WB
GTX646229	   OR2W3 antibody [HL5555]	Rec Rb mAb	Hu	WB
GTX637907	  ORP150 antibody [HL2006]	Rec Rb mAb	Hu, Ms, Rat	WB
GTX644843	   OXER1 antibody [HL4831]	Rec Rb mAb	Hu	WB
GTX647572	   P Glycoprotein antibody [GT391]	Rec Ms mAb	Hu, Ms	WB, IHC-P
GTX641111	   p130Cas (phospho Tyr165) antibody [HL3339]	Rec Rb mAb	Hu, Ms	WB
GTX641238	   p130Cas (phospho Tyr410) antibody [HL3401]	Rec Rb mAb	Hu, Ms	WB
GTX638869	    p27 Kip1 antibody [HL2508]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX645722	   P2Y4 antibody [HL5284]	Rec Rb mAb	Hu	WB
GTX639570	    p38 MAPK (phospho Thr180 / Tyr182) antibody [HL2741]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX639569	   p38 MAPK (phospho Tyr182) antibody [HL2740]	Rec Rb mAb	Hu	WB, IHC-P
GTX635797	    p38 MAPK antibody [HL1006]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	WB
GTX648087	   p40 antibody [HL6013]	Rec Rb mAb	Hu	WB, ICC/IF
GTX641465	   p53 (phospho Ser33) antibody [HL3559]	Rec Rb mAb	Hu	WB
GTX641560	   p53 (phospho Ser46) antibody [HL3589]	Rec Rb mAb	Hu	WB
GTX640490	   p53 (phospho Ser392) antibody [HL3047]	Rec Rb mAb	Hu	WB
GTX636395	    p53 (R175H Mutant) antibody [HL1129]	Rec Rb mAb	Hu	WB, ICC/IF
GTX638291	     p53 antibody [HL2242]	Rec Rb mAb	Cat, Dog, Hu, Ms	WB, ICC/IF, IHC-P
GTX636840	   p53R2 antibody [HL1400]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX635841	    p63 antibody [HL1027]	Rec Rb mAb	Hu, Ms, Zfsh	WB, ICC/IF, IHC-P, IHC-Wm
GTX641928	   p70 S6K (phospho Thr389) antibody [HL3738]	Rec Rb mAb	Hu	WB
GTX635621	      p70 S6K (phospho Thr421/Ser424) antibody [HL129]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX647218	   p73 antibody [HL5894]	Rec Rb mAb	Hu	WB, ICC/IF
GTX642041	   p90 RSK (phospho Ser380) antibody [HL3817]	Rec Rb mAb	Hu	WB, ICC/IF
GTX642333	   p90 RSK (phospho Thr573) antibody [HL4002]	Rec Rb mAb	Hu	WB
GTX642042	   p90 RSK1 (phospho Ser380) + RSK2 (phospho Ser386) antibody [HL3818]	Rec Rb mAb	Hu, Ms	WB
GTX641271	   p90 RSK1 (phospho Thr359) antibody [HL3429]	Rec Rb mAb	Hu	WB
GTX643556	    PAD4 antibody [HL4309]	Rec Rb mAb	Hu, Ms	WB, ICC/IF
GTX640408	   PAF Receptor antibody [HL2998]	Rec Rb mAb	Dog, Hu	WB
GTX638222	   PAI2 antibody [HL2218]	Rec Rb mAb	Hu, Ms	WB, IHC-P
GTX638484	  pan Cytokeratin antibody [HL2321]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX642052	   PAR3 antibody [HL3828]	Rec Rb mAb	Hu	ICC/IF
GTX636804	      PARP antibody [HL1364]	Rec Rb mAb	Hu, Ms	WB, ICC/IF, IHC-P, FCM, IP
GTX636983	    PARP3 antibody [HL1509]	Rec Rb mAb	Hu, Ms	WB, IHC-P
GTX648108	   PARP4 antibody [HL6034]	Rec Rb mAb	Hu	WB
GTX640615	   Paxillin (phospho Tyr31) antibody [HL3129]	Rec Rb mAb	Hu, Ms	WB, ICC/IF
GTX641204	   Paxillin (phospho Tyr118) antibody [HL3384]	Rec Rb mAb	Hu, Ms	WB
GTX641559	   Paxillin (phospho Ser273) antibody [HL3588]	Rec Rb mAb	Hu, Ms	WB
GTX644830	   Paxillin antibody [HL4818]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF
GTX645429	    PCBP2 antibody [HL5177]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX642583	   PCDHGB2 antibody [GT156]	Rec Ms mAb	Ms, Rat	WB
GTX637858	   PCNA antibody [HL1965]	Rec Rb mAb	Cat, Dog, Hu, ... ¹	WB, IHC-P, IHC-Fr, IHC-Wm
GTX643643	    PD1 antibody [HL4386]	Rec Rb mAb	Hu	WB, IHC-P
GTX646871	    PDGF Receptor alpha (D842V Mutant) antibody [HL5800]	Rec Rb mAb	Hu	WB
GTX641117	   PDGF Receptor alpha (phospho Tyr849) + PDGF Receptor beta (phospho Tyr857) antibody [HL3345]	Rec Rb mAb	Hu, Ms	WB
GTX645094	   PDGF Receptor alpha antibody [HL5001]	Rec Rb mAb	Hu, Ms	WB, ICC/IF
GTX640887	   PDGF Receptor beta (phospho Tyr740) antibody [HL3244]	Rec Rb mAb	Ms	WB
GTX641205	    PDGF Receptor beta (phospho Tyr751) antibody [HL3385]	Rec Rb mAb	Hu, Ms	WB, ICC/IF

¹ Cat, Dog, Hu, Japanese Medaka, Ms, Rat, Zfsh




LC3B antibody [HL2798]
(GTX639686)



Caspase 3 antibody [HL1370]
(GTX636810)




NRF1 antibody [HL3308]
(GTX641070)




Bax antibody [HL236]
(GTX635715)


Cat. No.		Product Name	Host & Clonality	Reactivity	Applications
GTX640692		PDGF Receptor beta (phospho Tyr1009) antibody [HL3174]	Rec Rb mAb	Hu, Ms	WB, ICC/IF
GTX635975		PD-L1 antibody [HL1041]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P, FCM
GTX640609		PEDF antibody [HL3123]	Rec Rb mAb	Dog, Hu	WB
GTX638548		PELO antibody [HL2345]	Rec Rb mAb	Hu, Japanese Medaka, Ms, Rat, Zfsh	WB, IHC-P, IHC-Wm
GTX638084		Periostin antibody [HL2112]	Rec Rb mAb	Hu, Ms	WB
GTX648919		Peripherin antibody [HL6471]	Rec Rb mAb	Hu, Ms	WB
GTX642201		Peroxisedoxin 6 antibody [HL3951]	Rec Rb mAb	Hu, Ms, Rat, Zfsh	WB, IHC-P
GTX638208		PFKFB4 antibody [HL2204]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX641916		PHACTR1 antibody [HL3726]	Rec Rb mAb	Ms, Rat	WB, IHC-P
GTX637407		PIGR antibody [HL1757]	Rec Rb mAb	Hu	WB
GTX641554		PKA R2 (phospho Ser99) antibody [HL3583]	Rec Rb mAb	Ms	WB
GTX636681		PKC (pan) (zeta phospho Thr410) antibody [HL1279]	Rec Rb mAb	Dm, Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX641580		PKC alpha (phospho Thr497) antibody [HL3605]	Rec Rb mAb	Hu, Ms	WB
GTX648872		PKC alpha antibody [HL6427]	Rec Rb mAb	Hu, Ms, Rat	WB
GTX636461		PKC theta (phospho Thr219) antibody [HL1157]	Rec Rb mAb	Hu	WB, ICC/IF, ELISA, Sandwich ELISA
GTX641233		PKC theta (phospho Ser695) antibody [HL3396]	Rec Rb mAb	Hu	WB
GTX636912		PKR (phospho Thr446) antibody [HL1439]	Rec Rb mAb	Hu	WB, ICC/IF
GTX637979		Plasminogen antibody [HL2071]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat, Zfsh	WB, ICC/IF, IHC-P
GTX639346		PLGF antibody [HL2682]	Rec Rb mAb	Hu	ELISA, Sandwich ELISA
GTX641191		PLK4 antibody [HL3371]	Rec Rb mAb	Hu, Ms	WB, IHC-P
GTX647176		PLVAP antibody [HL5858]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF
GTX637951		PODXL antibody [HL2047]	Rec Rb mAb	Cat, Hu, Ms, Rat	WB, IHC-P
GTX637647		Pogz antibody [HL1898]	Rec Rb mAb	Hu	WB
GTX638128		POLE antibody [HL2150]	Rec Rb mAb	Hu, Ms	WB, IHC-P
GTX641310		POLR3A antibody [HL3437]	Rec Rb mAb	Hu, Japanese Medaka, Ms, Rat, Zfsh	WB, ICC/IF, IHC-P
GTX636919		PP2A alpha antibody [HL1447]	Rec Rb mAb	Cat, Dm, Dog, Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX637433		PPP2R1A antibody [HL1781]	Rec Rb mAb	Hu, Ms, Rat	WB
GTX640828		PRAS40 (phospho Ser183) antibody [HL3189]	Rec Rb mAb	Hu	WB
GTX640987		PRAS40 (phospho Thr246) antibody [HL3298]	Rec Rb mAb	Hu, Ms	WB
GTX638919		PRG2 antibody [HL2547]	Rec Rb mAb	Hu, Ms	WB, ICC/IF, IHC-P
GTX639327		PSAT1 antibody [HL2663]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX639621		PVRL4 antibody [HL2758]	Rec Rb mAb	Hu	WB, FCM
GTX642098		PYK2 (phospho Tyr402) antibody [HL3860]	Rec Rb mAb	Hu, Ms	WB
GTX645213		Rab11A antibody [HL5043]	Rec Rb mAb	Hu, Rat	WB
GTX648441		RAB1B antibody [HL6331]	Rec Rb mAb	Hu, Ms, Rat	WB
GTX645495		RAB3C antibody [HL5239]	Rec Rb mAb	Hu, Ms, Rat	WB
GTX635981		RAB6A antibody [HL1047]	Rec Rb mAb	Dm, Hu, Ms, Rat, Zfsh	WB, IHC-P
GTX646104		RAB7A antibody [HL5462]	Rec Rb mAb	Hu, Ms	WB, IHC-P
GTX647254		RAB7L1 antibody [HL5929]	Rec Rb mAb	Hu, Ms, Rat	WB
GTX644762		RAB9B antibody [HL4758]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF
GTX641704		Rad17 (phospho Ser645) antibody [HL3676]	Rec Rb mAb	Hu	WB, ICC/IF
GTX641258		Rad21 antibody [HL3416]	Rec Rb mAb	Hu, Ms	WB, ICC/IF, IHC-P
GTX644763		RAB27A antibody [HL4759]	Rec Rb mAb	Dog, Hu, Ms, Rat	WB
GTX645073		RAB27B antibody [HL4987]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P
GTX646277		RAB31 antibody [HL5595]	Rec Rb mAb	Hu	WB
GTX70230		Rad51 antibody [14B4]	Rec Ms mAb	Chk, Hu, Ms, Rat	WB, ICC/IF, IHC-P, IP, IHC, PLA
GTX645970		Radixin antibody [HL5402]	Rec Rb mAb	Hu, Ms	WB, IHC-P
GTX647172		RALBP1 antibody [HL5854]	Rec Rb mAb	Hu	WB, ICC/IF
GTX638564		RanBP1 antibody [HL2358]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX638871		RanBP16 antibody [HL2510]	Rec Rb mAb	Hu, Ms	WB, ICC/IF




PARP antibody [HL1364]
(GTX636804)     



NfκB p65 antibody [HL4229]
(GTX643062)    



Bcl-XL antibody [HL2039]
(GTX637940)    



SQSTM1 / P62 antibody [HL1113]
(GTX636328)    


Cat. No.	Product Name	Host & Clonality	Reactivity	Applications
GTX638637	  RAP1A antibody [HL2415]	Rec Rb mAb	Hu	WB
GTX636434	   RAS (G12A mutant) antibody [HL1153]	Rec Rb mAb	Hu	WB
GTX635362	    RAS (G12D Mutant) antibody [HL10]	Rec Rb mAb	Hu, Ms	WB, IHC-P
GTX635623	    RAS (G12V Mutant) antibody [HL169]	Rec Rb mAb	Hu	WB, IHC-P
GTX646283	  RAS (Q61H mutant) antibody [HL5601]	Rec Rb mAb	Hu	WB, IHC-P (cell pellet)
GTX640939	  Rb (phospho Ser780) antibody [HL3285]	Rec Rb mAb	Hu	WB
GTX638814	 Rb antibody [HL2465]	Rec Rb mAb	Hu	WB, ICC/IF
GTX648084	  RCC1 antibody [HL6010]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX638267	   RelB antibody [HL2222]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P, IP
GTX640885	 RFC1 antibody [HL3242]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF
GTX640397	 RGS2 antibody [HL2987]	Rec Rb mAb	Ms, Rat, Zfsh	WB, IHC-P
GTX645790	  RICTOR antibody [HL5325]	Rec Rb mAb	Hu, Ms, Rat	WB
GTX638299	  RIP antibody [HL2250]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P
GTX638546	 RNase H1 antibody [HL2343]	Rec Rb mAb	Dog, Hu, Ms	WB, IHC-P
GTX638543	 RNase L antibody [HL2340]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX638350	  RNF31 antibody [HL2295]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX644685	  RNF43 antibody [HL4684]	Rec Rb mAb	Hu, Ms	WB
GTX637125	  ROCK1 + ROCK2 antibody [HL1632]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX636103	   ROCK2 (phospho Ser1366) antibody [HL1071]	Rec Rb mAb	Hu, Ms	WB, ICC/IF
GTX648158	   ROCK2 antibody [HL6079]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF
GTX647207	  RPA32 (phospho Ser4) antibody [HL5883]	Rec Rb mAb	Hu	WB, ICC/IF
GTX641864	  RPA32 (phospho Thr21) antibody [HL3692]	Rec Rb mAb	Hu, Ms	WB, ICC/IF
GTX644832	 RPL3 antibody [HL4820]	Rec Rb mAb	Dog, Hu, Ms, Rat	WB, IHC-P
GTX638494	  RPL32 antibody [HL2331]	Rec Rb mAb	Dog, Hu, ... ^①	WB, ICC/IF, IHC-P, IHC-Wm
GTX641110	 RPS2 antibody [HL3338]	Rec Rb mAb	Hu, Ms, Rat	WB
GTX642012	  RPS6 (phospho Ser235) antibody [HL3793]	Rec Rb mAb	Hu, Ms	WB
GTX642048	  RPS6 (phospho Ser240/Ser244) antibody [HL3824]	Rec Rb mAb	Hu, Ms	WB, ICC/IF
GTX636579	 RRM1 antibody [HL1231]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX636634	  RRM2 antibody [HL1235]	Rec Rb mAb	Hu, Ms, Zfsh	WB, ICC/IF, IHC-P, IHC-Wm
GTX641269	   RSK (pan) (phospho Ser221 / Ser227 / Ser218 / Ser232) antibody [HL3427]	Rec Rb mAb	Hu, Ms	WB
GTX640988	   RSK3 (phospho Thr356/Ser360) antibody [HL3299]	Rec Rb mAb	Hu, Ms	WB, ICC/IF
GTX648350	  S100A4 antibody [HL6258]	Rec Rb mAb	Hu, Ms	WB, IHC-P
GTX648148	  S100A9 antibody [HL6069]	Rec Rb mAb	Hu, Ms	WB
GTX645468	  SA1 antibody [HL5212]	Rec Rb mAb	Hu	WB, ICC/IF
GTX636097	    SARM1 antibody [HL1066]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX636098	   SDHA antibody [HL1067]	Rec Rb mAb	Dm, Hu, ... ^②	WB, ICC/IF
GTX647598	  SEC61B antibody [HL5984]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX637564	   SEK1 / MKK4 antibody [HL1839]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX646333	  SHARP2 antibody [HL5644]	Rec Rb mAb	Hu	IHC-P
GTX645456	  SHIP1 antibody [HL5204]	Rec Rb mAb	Hu, Ms	WB, ICC/IF
GTX637909	  SIGIRR antibody [HL2008]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX640659	  SIRP alpha antibody [HL3149]	Rec Rb mAb	Hu	WB
GTX642513	 SLBP antibody [HL4074]	Rec Rb mAb	Hu	WB, ICC/IF
GTX637889	  SLC1A5 antibody [HL1988]	Rec Rb mAb	Cat, Hu, Ms	WB, ICC/IF, IHC-P
GTX638277	  SLC27A2 antibody [HL2232]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX648130	  SLC39A4 antibody [HL6051]	Rec Rb mAb	Hu	WB
GTX639352	  SLC5A8 antibody [HL2687]	Rec Rb mAb	Ms	WB, IHC-P
GTX641461	  SMAD2 (phospho Ser245/Ser250/Ser255) antibody [HL3555]	Rec Rb mAb	Hu	WB
GTX641199	  SMAD2 (phospho Ser465 / Ser467) + SMAD3 (phospho Ser423 / Ser425) antibody [HL3379]	Rec Rb mAb	Hu, Ms	WB

^① Dog, Hu, Japanese Medaka, Ms, Rat, Zfsh


^② Dm, Hu, Japanese Medaka, Ms, Rat, Zfsh




Ki67 antibody [HL4117]
(GTX642563)



PCNA antibody [HL1965]
(GTX637858)




MCM7 antibody [HL4072]
(GTX642511)




c-Myc antibody [HL6095]
(GTX648180)


Cat. No.	Product Name	Host & Clonality	Reactivity	Applications
GTX643573	SMAD2 / SMAD3 (phospho Thr8) antibody [HL4326]	Rec Rb mAb	Hu, Ms	WB
GTX636329	SMARCA2 antibody [HL1114]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX635798	SMC8 antibody [HL1007]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX638370	SNA1 antibody [HL2303]	Rec Rb mAb	Hu, Ms	WB, ICC/IF
GTX638270	SNRK antibody [HL2225]	Rec Rb mAb	Hu, Rat	WB
GTX636934	SOD2 (Acetyl Lys68) antibody [HL1460]	Rec Rb mAb	Hu, Ms	WB, Dot
GTX636957	SOD2 antibody [HL1483]	Rec Rb mAb	Hu, Japanese Medaka, Zfsh	WB, IHC-P
GTX635461	Sodium/Potassium ATPase alpha 1 antibody [HL114]	Rec Rb mAb	Dm, Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX639345	Somatostatin receptor 3 antibody [HL2681]	Rec Rb mAb	Hu, Ms	WB, IHC-P
GTX636505	SOX2 antibody [HL1193]	Rec Rb mAb	Hu, Ms, Rat, Cat	WB, ICC/IF, IHC-P, IP
GTX645986	SPINK4 antibody [HL5418]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX641562	SQSTM1 / P62 (phospho Ser272) antibody [HL3591]	Rec Rb mAb	Hu, Ms	WB
GTX645296	SQSTM1 / P62 (phospho Ser349) antibody [HL5102]	Rec Rb mAb	Hu, Rat	WB
GTX636328	SQSTM1 / P62 antibody [HL1113]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX641636	Src (phospho Tyr530) antibody [HL3645]	Rec Rb mAb	Hu	WB
GTX647232	SSB antibody [HL5907]	Rec Rb mAb	Hu, Ms	WB, IHC-P
GTX645293	ST2 antibody [HL5099]	Rec Rb mAb	Hu	WB
GTX641342	ST3GAL2 antibody [HL3457]	Rec Rb mAb	Hu	WB
GTX638481	ST6GAL1 antibody [HL2318]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX642049	STAT1 (phospho Tyr701) antibody [HL3825]	Rec Rb mAb	Hu	ICC/IF
GTX643031	STAT2 (phospho Tyr690) antibody [HL4208]	Rec Rb mAb	Hu, Ms, Rat	WB
GTX642311	STAT2 antibody [HL3985]	Rec Rb mAb	Hu	WB
GTX636938	STAT3 (phospho Tyr705) antibody [HL1464]	Rec Rb mAb	Hu, Rat	WB
GTX636400	STAT3 antibody [HL1134]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	WB, IHC-P
GTX641431	STX17 antibody [HL3534]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX637784	SUN1 antibody [HL1947]	Rec Rb mAb	Hu, Ms	WB, ICC/IF, IHC-P
GTX641194	Suppressor of Fused antibody [HL3374]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX641550	SYT13 antibody [HL3579]	Rec Rb mAb	Hu	WB
GTX637311	TAB2 antibody [HL1712]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX647224	TACSTD2 antibody [HL5900]	Rec Rb mAb	Hu	WB, IHC-P
GTX637302	TAL1 antibody [HL1703]	Rec Rb mAb	Hu	WB
GTX645108	TAS2R3 antibody [HL5015]	Rec Rb mAb	Hu	WB
GTX645491	TAS2R7 antibody [HL5235]	Rec Rb mAb	Hu	WB
GTX644970	TAS2R40 antibody [HL4905]	Rec Rb mAb	Hu	WB, ICC/IF
GTX644971	TAS2R50 antibody [HL4906]	Rec Rb mAb	Hu	WB, ICC/IF, FCM
GTX637791	TCF3 / E2A antibody [HL1954]	Rec Rb mAb	Hu	WB, ICC/IF
GTX645771	TFEB (phospho Ser211) antibody [HL5306]	Rec Rb mAb	Hu	ICC/IF
GTX636089	TFEB antibody [HL1065]	Rec Rb mAb	Hu, Ms	WB, IHC-P
GTX640355	TFF1 antibody [HL2968]	Rec Rb mAb	Hu	WB
GTX640232	TGF beta Receptor III antibody [HL2895]	Rec Rb mAb	Hu	WB, ICC/IF
GTX636411	Thrombospondin 2 antibody [HL1140]	Rec Rb mAb	Hu	WB, ICC/IF, ELISA, Sandwich ELISA
GTX636635	Thymidylate synthase antibody [HL1236]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX639116	TIGIT antibody [HL2657]	Rec Rb mAb	Hu	WB, IHC-P
GTX638477	TLR3 antibody [HL2314]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	WB, IHC-P
GTX637058	TNF alpha antibody [HL1579]	Rec Rb mAb	Hu, Ms	WB, ELISA
GTX637401	TNK1 antibody [HL1751]	Rec Rb mAb	Hu	WB
GTX637534	TOX antibody [HL1809]	Rec Rb mAb	Hu	WB
GTX642127	TRAF6 antibody [HL3886]	Rec Rb mAb	Hu	ICC/IF
GTX636672	Transgelin antibody [HL1270]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX642033	TRAP1 antibody [HL3809]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF




Histone H2A.X antibody [HL3515]
(GTX641412)



TACSTD2 antibody [HL5900]
(GTX647224)



PDGF Receptor beta
(phospho Tyr751) antibody
[HL3385] (GTX641205)



xCT / SLC7A11 antibody [HL5661]
(GTX646470)

Cat. No.	Product Name	Host & Clonality	Reactivity	Applications
GTX635872	TREM2 (R47H Mutant) antibody [HL1032]	Rec Rb mAb	Hu	WB
GTX637386	TREM2 antibody [HL1738]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P, IP, ELISA, Sandwich ELISA
GTX648964	TRP2 antibody [HL6508]	Rec Rb mAb	Hu	WB
GTX638190	Tryptophanyl tRNA synthetase antibody [HL2186]	Rec Rb mAb	Hu, Ms, Rat, Zfsh	WB, ICC/IF, IHC-P
GTX640361	TSC2 antibody [HL2974]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX70255	TSG101 antibody [4A10]	Rec Ms mAb	Dog, Hm, Hu, Mk, Ms, Pig, Rat, Zfsh	WB, ICC/IF, IHC-P, FCM, IP, ELISA, EM, IHC, ... ¹
GTX644761	TXA2R antibody [HL4757]	Rec Rb mAb	Hu	WB
GTX640840	TXNDC5 antibody [HL3200]	Rec Rb mAb	Hu	WB, IHC-P
GTX638210	ULK2 antibody [HL2206]	Rec Rb mAb	Hu, Rat	WB, ICC/IF
GTX646215	UPK1A antibody [HL5541]	Rec Rb mAb	Ms	WB, IHC-P
GTX644758	UPK1B antibody [HL4754]	Rec Rb mAb	Hu, Ms	WB
GTX647577	USP11 antibody [HL5963]	Rec Rb mAb	Hu, Rat	WB, ICC/IF
GTX646149	USP7 antibody [HL5485]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P
GTX648072	VCAM1 / CD106 antibody [HL5998]	Rec Rb mAb	Hu	WB, IHC-P
GTX641193	VDAC1 antibody [HL3373]	Rec Rb mAb	Hu, Ms	WB
GTX648354	VE-Cadherin antibody [HL6262]	Rec Rb mAb	Hu, Ms	WB
GTX642018	VEGF Receptor 1 antibody [HL3799]	Rec Rb mAb	Hu	IHC-P
GTX645284	VEGF Receptor 3 antibody [HL5094]	Rec Rb mAb	Hu	WB, ICC/IF
GTX637405	VEGFA antibody [HL1755]	Rec Rb mAb	Cat, Dog, Hu, Ms	WB, ICC/IF
GTX638344	VEGFC antibody [HL2289]	Rec Rb mAb	Hu	WB, ICC/IF
GTX640850	VGLuT1 antibody [HL3210]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P, IHC-Fr
GTX637555	Villin antibody [HL1830]	Rec Rb mAb	Hu, Ms	WB, IHC-P
GTX636980	Vimentin antibody [HL1506]	Rec Rb mAb	Cat, Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX637857	Vinculin antibody [HL1964]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX639325	Vitamin D Receptor antibody [HL2661]	Rec Rb mAb	Hu	WB
GTX643590	VPS4B antibody [HL4340]	Rec Rb mAb	Hu	WB
GTX635821	VPS35 antibody [HL1017]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat, Zfsh	WB, ICC/IF, IHC-P, IP
GTX646070	VTA1 antibody [HL5428]	Rec Rb mAb	Hu	WB, ICC/IF
GTX641318	WAVE1 antibody [HL3445]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P, IHC-Fr
GTX637890	WDFY1 + WDFY2 antibody [HL1989]	Rec Rb mAb	Hu	WB
GTX644984	WEE1 antibody [HL4911]	Rec Rb mAb	Hu, Ms	WB
GTX637638	Wnt4 antibody [HL1889]	Rec Rb mAb	Hu, Ms, Zfsh	WB, ICC/IF
GTX637377	WRN antibody [HL1729]	Rec Rb mAb	Hu	WB, IHC-P
GTX646470	xCT / SLC7A11 antibody [HL5661]	Rec Rb mAb	Hu	WB, ICC/IF
GTX640231	XPC antibody [HL2894]	Rec Rb mAb	Hu, Ms, Rat	WB, IHC-P
GTX641911	YAP1 (phospho Ser61) antibody [HL3721]	Rec Rb mAb	Hu	WB
GTX637643	YAP1 antibody [HL1894]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P
GTX637587	YTHDC2 antibody [HL1862]	Rec Rb mAb	Hu, Ms	WB, ICC/IF, IHC-P
GTX643021	YY1 antibody [HL4198]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX636914	ZC3H12A antibody [HL1442]	Rec Rb mAb	Cat, Hu	WB
GTX637633	ZC3HAV1 antibody [HL1884]	Rec Rb mAb	Hu	WB, ICC/IF, IHC-P
GTX638294	ZEB1 antibody [HL2245]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
GTX641379	ZER1 antibody [HL3491]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF
GTX636491	ZO-1 antibody [HL1185]	Rec Rb mAb	Cat, Dog, Hu, Ms, Rat	WB, ICC/IF
GTX648337	ZO-2 antibody [HL6245]	Rec Rb mAb	Hu, Ms, Rat	WB

¹ WB, ICC/IF, IHC-P, FCM, IP, ELISA, EM, IHC, IHC (Free Floating)