


The Pre-Metastatic Niche



Recombinant | KO/KO Validation | Orthogonal Validation | Protein Overexpression | HiStoMAX

ECM Remodeling

 Collagen II antibody [HL1772] (GTX637424)	 Collagen III antibody [HL1906] (GTX637655)	 Collagen I + Collagen II + Collagen III antibody [HL2048 + HL1907] (GTX638633)	 COL4A1 antibody [HMV327] HistoMAX™ (GTX640559)	 Elastin antibody [MSVA-648R] HistoMAX™ (GTX04466)	 FAK (phospho Tyr576) antibody [HL127] (GTX635714)	 S100A12 antibody [MSVA-812M] HistoMAX™ (GTX04449)	 S100P antibody [MSVA-480R] HistoMAX™ (GTX04398)	 S100 beta antibody [MSVA-490R] HistoMAX™ (GTX04404)
---	--	--	--	---	---	---	---	---

Inflammation & Immunosuppression & Pre-Metastatic Niche Formation

 TNF alpha antibody [HL2361] (GTX638567) TNF alpha antibody [HL2363] (GTX638569)	 VEGFC antibody [HL2289] (GTX638344)	References: <ol style="list-style-type: none"> 1. Annu Rev Pathol. 2021 Jan 24;16:409-432. 2. Nat Rev Cancer. 2017 May;17(5):302-317. 3. Nat Rev Immunol. 2023 Aug;23(8):522-538. 4. Cancer Cell. 2023 Mar 13;41(3):546-572.
--	---	---

