


Organelle Markers

Organelle Markers


Nucleus

HDAC2 antibody [HL1545] (GTX637020)




Nucleolus

PAF49 antibody (GTX102175)




Microtubules

alpha Tubulin antibody (GTX112141)




Exosome

CD63 antibody (GTX135220)




Nuclear Envelope

Lamin A + C antibody (GTX101127)




Endoplasmic Reticulum

GRP94 antibody (GTX103232)




Actin

beta Actin antibody (GTX109639)




Mitochondria

COX4 antibody (GTX114330)




Centromere

gamma Tubulin antibody [HL1175] (GTX636480)



Cytoplasm

Hsp70 antibody (GTX111088)




Antibody List

Organelle	Cat. No.	Product Name	Clonality	Reactivity	Applications
Actin	GTX109639	beta Actin antibody	Rb pAb	Hu, Ms, Rat, Zfish, Dog, Cat, Dog, Hm, Chk, Pig, Mf, C. elegans, C. albicans, E. coli, Hrs, Rice, Xenopus, Drosophila, Plant, Worm	WB, ICC/IF, IHC-P, IHC-Fr, IP
Cytoplasm	GTX111088	Hsp70 antibody	Rb pAb	Hu, Ms, Rat, Mf, C. elegans	WB, ICC/IF, IHC-P, IP
Plasma Membrane	GTX132646	pan Cadherin antibody	Rb pAb	Hu, Ms	WB, ICC/IF, FCM
	GTX635461	Sodium/Potassium ATPase alpha 1 antibody [HL114]	Rec Rb mAb	Hu, Ms, Rat, Dm	WB, ICC/IF, IHC-P
	GTX103232	GRP94 antibody	Rb pAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P, IP
Endoplasmic Reticulum	GTX637077	Calnexin antibody [HL1598]	Rec Rb mAb	Hu, Ms, Rat, Cat, Dog	WB, IHC-P
	GTX16154	GOLPH2 antibody	Rb pAb	Hu	WB, ICC/IF, IHC-P
Golgi Apparatus	GTX104255	GOLGA5 antibody [N2C2], Internal	Rb pAb	Hu, Ms, Hm	WB, ICC/IF, IHC-P
Junction (Tight Junction)	GTX636491	ZO-1 antibody [HL1185]	Rec Rb mAb	Hu, Ms, Rat, Cat, Dog	WB, ICC/IF
Junction (Gap Junction)	GTX133022	Connexin 43 antibody	Rb pAb	Hu, Ms	WB, ICC/IF, IHC-P
Centrosome	GTX636480	gamma Tubulin antibody [HL1175]	Rec Rb mAb	Hu, Ms, Rat, Zfish, Rb, Dm, Cat, Dog, Hm, AGMR	WB, ICC/IF, IHC-Fr
	GTX107155	ZWINT antibody	Rb pAb	Hu	WB, ICC/IF
Kinetochore (M Phase)	GTX02068	Hec1 antibody [9G3.23]	Ms mAb	Hu, Ms, Hm, Kg	WB, ICC/IF, IHC-Fr, Fcm, IP, PLA
	GTX634336	LAMP1 antibody [GT25212]	Ms mAb	Hu	WB, ICC/IF, IHC-Fr
Lysosome	GTX103214	LAMP2 antibody	Rb pAb	Hu	WB, ICC/IF, IHC-P
	GTX634169	EEA1 antibody [GT10811]	Ms mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
Endosome	GTX132548	RAB7A antibody	Rb pAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
	GTX119095	RAB11B antibody [N1C3]	Rb pAb	Hu, Ms, Rat	WB, ICC/IF, IHC-Fr
Exosome	GTX644095	CD9 antibody [HL4525]	Rec Rb mAb	Hu, Rat	WB, IHC-P, ELISA, Sandwich ELISA
	GTX643611	CD63 antibody [HL4361]	Rec Rb mAb	Hu, Ms, Cat	WB, IHC-P
	GTX101766	CD81 antibody	Rb pAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P, FCM, Blocking
Autophagosome	GTX02055	TSGL1 antibody [4A10]	Ms mAb	Hu, Ms, Rat, Zfish, Dog, Hm, Pig, Mf	WB, ICC/IF, IHC-P, FCM, IP, ELISA, EM, IHC, IHC (Free Floating)
	GTX127375	LC3B antibody	Rb pAb	Hu, Ms, Rat, Guinea pig, Pig, Mosq	WB, ICC/IF, IHC-P, IHC-Fr, FCM, IP, EM
	GTX636328	SQSTM1 / P62 antibody [HL1113]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
Microtubule / Spindle	GTX628802	alpha Tubulin antibody [GT114]	Ms mAb	Hu, Ms, Rat, Zfish, Goat, Dm, Crucian carp	WB, ICC/IF, IHC-P, IHC-Fr, FCM
Microtubule (Inner Membrane)	GTX114330	COX4 antibody	Rb pAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P, IHC-Fr
Microtubule (Outer Membrane)	GTX101105	Prohibitin antibody	Rb pAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P, IP
	GTX636967	TOMM20 antibody [HL1493]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
Nuclear Envelope	GTX631209	FS1 antibody [GT9810]	Ms mAb	Hu, Ms, Rat	WB, ICC/IF
	GTX101127	Lamin A + C antibody	Rb pAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P, IP
	GTX102128	POM121 antibody [N2N3]	Rb pAb	Hu, Ms	WB, ICC/IF, IHC-P
Nucleolus	GTX102175	PAF49 antibody	Rb pAb	Hu	WB, ICC/IF, IHC-P
	GTX108749	RPA70 antibody [C1C3]	Rb pAb	Hu, Ms	WB, ICC/IF, IHC-P, IP
Nucleus	GTX100539	PCNA antibody	Rb pAb	Hu, Ms, Rat, Dm, Golden Syrian Hamster	WB, ICC/IF, IHC-P, IP
	GTX100513	HDAC1 antibody	Rb pAb	Hu, Ms, Rat, Zfish	WB, ICC/IF, IHC-P, IHC-Fr, IHC-Wm, IP, CHIP assay
	GTX637020	HDAC2 antibody [HL1545]	Rec Rb mAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P
	GTX109749	Vinculin antibody [N3C1], Internal	Rb pAb	Hu, Ms, Rat, Pig	WB, IHC-P
	GTX102156	Cofilin 1 antibody	Rb pAb	Hu, Ms, Rat	WB, ICC/IF, IHC-P, IHC-Fr
	GTX104698	Cyclophilin A antibody [C1C3]	Rb pAb	Hu, Ms, Rat, Hm	WB, ICC/IF, IHC-P
Internal Controls	GTX122148	Histone H3 antibody	Rb pAb	Hu, Ms, Rat, Dm, Pig, Mf, Fungi, Golden Syrian Hamster, Rice	WB, ICC/IF, IHC-P, IP, CHIP assay, IHC

More Antibodies


Tubulin

alpha Tubulin antibody [GT114] (GTX628802)




Actin

alpha Skeletal Muscle Actin antibody (GTX101362)




Mitochondria

COX4 antibody (GTX101499)




Nucleus

Prohibitin antibody (GTX101105)




Internal Controls

Histone H3 antibody (GTX122148)




Nucleolus

NOL5A antibody (GTX132073)




Nuclear Envelope

POM121 antibody [N2N3] (GTX102128)




Cytoplasm

GAPDH antibody (GTX100118)




Endoplasmic Reticulum

Calnexin antibody [C3], C-term (GTX108660)




Golgi Apparatus

GOLPH2 antibody (GTX118154)




Lysosome

LAMP2 antibody (GTX103214)



Autophagosome

LC3B antibody (GTX127375)



International
Tel : 866.3.6208988 Fax : 866.3.6208989
E-mail : international@genetex.com

USA
Tel : 1.949.553.1900 Fax : 1.949.309.2888
E-mail : info@genetex.com

Europe
Tel : 49.8161.2429021 Fax : 49.8161.2429069
E-mail : europa@genetex.com

